

Product datasheet for **RR207349L4V**

Phkg1 (NM_031573) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Phkg1 (NM_031573) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Phkg1
Synonyms:	Phkg; PHKIN01
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_031573
ORF Size:	1164 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR207349).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_031573.1 , NP_113761.1
RefSeq Size:	1388 bp
RefSeq ORF:	1167 bp
Locus ID:	29353
UniProt ID:	P13286
Cytogenetics:	12q13
Gene Summary:	catalytic subunit of phosphorylase kinase complex; plays a role in glycogen metabolism [RGD, Feb 2006]



[View online »](#)