

## Product datasheet for **RR207357**

### Ldha (NM\_017025) Rat Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Ldha (NM_017025) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Ldha
Synonyms:	CDK1; Ldh1
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin
ORF Nucleotide Sequence:	>NM_017025 ORF sequence, RR207357 may differ due to SNPs. Blue=ORF Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGGCAGCCCTCAAGGACCAGCTGATTGTGAATCTTCTTAAGGAAGAACAGGTCCCCAGAACAAGATT
ACAGTTGTTGGGGTTGGTGTGCTGTTGGCATGGCTTGTGCCATCAGTATCTTAATGAAGGACTTGGCTGAT
GAGCTTGGCCCTTGTGATGTCATAGAAGATAAGCTAAAGGGAGAGATGATGGATCTTCAGCATGGCAGC
CTTTTCTTAAGACACCAAAAATTGTCTCCAGCAAAGATTATAGTGTGACTGCAAACCTCCAAGCTGGTC
ATTATCACCGCGGGGCCCGTCAGCAAGAGGGAGAGAGCCGGCTCAATTTGGTCCAGCGAAACGTGAAC
ATCTTCAAGTTTCATCATTCCAAATGTTGTGAAATACAGTCCACAGTGCAAACCTGCTCATCGTCTCAAAC
CCAGTGGATATCTTGACCTACGTGGCTTGGGAAGATCAGCGGCTTCCCCAAAAACAGAGTTATTGGAAGT
GGTTGCAATCTGGATTCCGGCTCGGTTCCGTTACCTGATGGGAGAAAGGCTGGGAGTTCACTCACTGAGC
TGTCACGGGTGGGTCTGGGAGAGCATGGCGACTCCAGTGTGCCTGTGTGGAGTGGTGTGAATGTCGCC
GGCGTCTCCCTGAAGTCTCTGAACCCGAGCTGGGCACGGATGCAGACAAGGAGCAGTGGAAAGGATGTG
CACAAGCAGGTGGTTGACAGTGCATACGAAGTATCAAGCTGAAAGGTTACACATCCTGGGCCATTGGC
CTCTCCGTGGCAGACTTGGCCGAGAGCATAATGAAGAACCTTAGGCGGGTGCATCCCATTCCACCATG
ATTAAGGGTCTCTATGGAATCAAGGAGGATGTCTTCTCAGCGTCCCATGTATCCTGGGACAAAATGGA
ATCTCAGATGTTGTGAAGGTGACACTGACTCCTGACGAGGAGGCCCGCCTGAAGAAGAGTGCAGATACC
CTCTGGGGAATCCAGAAGGAGCTGCAGTTC
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```




[View online »](#)

**Protein Sequence:** >Peptide sequence encoded by RR207357  
 Blue=ORF Red=Cloning site Green=Tag(s)


MAALKDQLIVNLLKEEQVPQNKITVVGVGAVGMACAISILMKDLADELALVDVIEDKLGEMMDLQHGS  
 LFLKTPKIVSSKDYSVTANSKLVITAGARQQEGESRLNLVQRNVNIFKFIIPNVVKYSPQCKLLIVSN  
 PVDILTYYAWKISGFPKNRVIGSGCNLDSARFRYLMGERLGVHPLSCHGWVLGEHGDSSVPVWSGVNVA  
 GVSLKSLNPQLGTDADKEQWKDVHKQVVD SAYEVIKLGKGYT SWAIGL SVADLAESIMKNLRRVHP ISTM  
 IKGLYGIKEDVFLSVPCILGQNGISDVVKVLTLPDEEARLKKSADTLWGIQKELQF  
 TRTRPLEQKLISEEDLAANDILDYKDDDDKV

**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**



**Plasmid Map:**



**ACCN:** NM\_017025  
**ORF Size:** 996 bp

<b>OTI Disclaimer:</b>	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq Size:</b>	1598 bp
<b>RefSeq ORF:</b>	999 bp
<b>Locus ID:</b>	24533
<b>UniProt ID:</b>	<a href="#">P04642</a>
<b>Cytogenetics:</b>	1q22
<b>MW:</b>	36.5 kDa
<b>Gene Summary:</b>	mRNA expression in fibroblasts increases in response to induction by epidermal growth factor or serum [RGD, Feb 2006]