

Product datasheet for **RR207389**

Muc20 (NM_001107089) Rat Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Muc20 (NM_001107089) Rat Tagged ORF Clone
Tag: Myc-DDK
Symbol: Muc20
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RR207389 representing NM_001107089
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGATGAGCTCTGTCTGGGGCTGGCTGTACCGCTCCTTGCTTCTACTGGAAGTTGGAGTCAGTGTGA
GCTCTCCAGGGACATTGGCACCAAGTGACAGCCCCAGAGCTGGGATGACCACTGTCAAAGCTGGCACAGT
GAGTGACCCCATAGAAGCCGCTCTTGGCACCCCTTGTACTGATGACAGCTCTGAAGAGACAAGGAAGATC
ACAATTGACCTCCTGACTTTGGCACACCCCTCCATGGGAGTGGAGCAGCTGTCTTCAGAGAGCAGCTCCT
CCTCTGACAGCTCAGCTGGGGTCTCAGCAGCTCTCAGGCCCTGGGACCCGACAGTACAACCTCCAGCCAA
GGACTCAGTTGCCTTCAGCATCTCCCATATTAAGCTCACCACTGCATCACGAAATAGAAAACAACCATC
ACCATCTCTGGGACCCAGACGCGAGCCGAGCCCTGCAGAAGCAATGACAGCTTTATCTGCTTCTGAGA
TGTTAACTTTGCCTCCACTGAGGCAAAGCCAGTTGTCCCTAAGACCACAAGCTCATCTGGGATCCT
GTCAACTGCCACTACCCTGGCCCTTGCACCACCCTTGAGGGTACTCCCTGCTAGTGAATCACAGAA
AGCGAAATAGCAGTGGCCAGACTCCTACCTCCAGTGGAACTTCGGCTACAGGAGGAACCCAGGTGACAG
TTAGAAGGAACCCACTGGAAGACTTTCAGCACTCTCTAGTGAGACACAAAGCCACTGAGTTTTTTGG
TACCATCACAGTCCCCACGGTTGCTGGATCAACAGTAGGAGAAGCAACCTCTCTGTGAGTCCACAGCT
CTAGACAGCAGTCTCTCAGCAGTGGTCACCACCAAGGGCTCTGCTTCTCAGAGACTTTAACACAGACA
ATACAACCAACAGCTCCTTCTCACAGGGAGTCGCCCTCCTTCTTAATCTATTGACAACAGCTAGCAC
CAGTGAAGAACAACAGTCACTTTAACCAAGACCACAGCCTCACCGAAGACCCCGATGAATCCCGCACCC
ACTGCTTGGACAAGGAAAACACAGAGCACGATCCAGGTATAAATGGAGGCTTCTCCTTGTACGGCTGA
CCGTGGCTCCCTAAGGACCTCACTGAGCACGTACCAGAGAGAAGCTGATGCACCAGCTCCGCCGTGA
ACTGCACACCCATATGCCACTTGTCCAAGTGTCTTACTGAGCATCAGGAGGAGC

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)


Protein Sequence: >RR207389 representing NM_001107089
 Red=Cloning site Green=Tags(s)

MMSSVWGLAVPLLVFYWKVGVSVSSPGTLAPSDSPRAGMTTVKAGTVSDPIEAVLGLTCTDSSSEETRKI
 TIDLLLAHPMSMGVEQLSSESSSSDSSAGVLSQQALGPDSTTPAKDSVAFSISHIKLTCITEIETTI
 TISGTPDASRSPAEAMTALSAEMLTLPSTEAKPVVPKTTSSSGILSTATTLALATTLEGLPASGITE
 SEIAVAQTPTSSGTSATGGTQVTVRRNPLEDTSALSSETQSHTEVFGTITVPTVAGSTVGEATSLVSSTA
 LDSSLSAVVTTKGSASSETLTDNTTNSFLTGSRPPSLIYSTTASTSERTNVTLTKTASPKTPMNPAP
 TAWTRKTEHDPGINGGFLVRLTVASPKDLTEHVTREKLMHQLRRELHTHMPLVQVSLLSIRRS


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001107089

ORF Size:	1245 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001107089.2 , NP_001100559.1
RefSeq Size:	1576 bp
RefSeq ORF:	1248 bp
Locus ID:	303886
Cytogenetics:	11q22
MW:	42.9 kDa