

Product datasheet for **RR207396**

Gpr108 (NM_199399) Rat Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Gpr108 (NM_199399) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Gpr108
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RR207396 representing NM_199399
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGGCAGTGAGCGAGAGGAGGGGGCTCAGCGCGAGAGCCAGCGCAGTGTGATGGAATACTTGTGCG
 TGCTGGTGCTGATGCTGAGCGGCTGCTCAGGCCGTATCCACCGGTGACTCTCACAGGGGAGAAGCGAGC
 TGACATCCAACCTGAACAGCTTTGGATTCTACACCAATGGTTCCTGGAAGTAGAGCTGAGCTTACTGCGA
 CTCAGCCTCCAGGAGACAGAAGACAAGTTTCAAAGGTGGGATTGAGCTTGGCCGAGTTCGGTCCGGCA
 GTGTTGCTCCTACTCGAGCCGAATCCCATGAATGCCCTCTGGAGAGAAACAGCAGCAACTTCTGGT
 CCTGTTCTCATCAATATCAAGGACCTGCAGGTCCAAGTCCGAAAATATGGGGAACAGAAGTTGTTTCATC
 TCCCTGGCTCCTCCCTGAAGTCCCTCCAGTCAGGGCCCCGAGCCAGACCAACAGGCACCCCCA
 AGGACAACCATGTCATTACCCATCTCCGAAAAGATGTCTGCAGTCAAGGAAGATCAGGCAAAGCTTAC
 AGTGCCCAAGTTTCCGGGGATAAAGTCTCCAGCGGGACACAGGCACTCATCTGATGGCCAGCCACAA
 AGCCAGCCACCCACAGGGGGCCAGTGGGAAGGAGAAGGATCTGGTTTTGGGACTGGGCCACCTCAATG
 ACTCTACAACCTCAGTTTCCACATTGTGATCGGCTCCAGGGCTGAGGAAGGCCAGTACAGCCTCAACT
 CCACAACCTGTTACAACACAATACCAGGCCAGGAGCAGCGTTCGACCTCACTGTAATGATCCGAGAGAAG
 AATCCAGAAGGCTTCTGTGAGCAGTGAATCCCTCTTCAAACCTCTATCTAATCATGTCCGCTGCT
 TCTTGGCCGCCGGCATCTTCTGGGTGTCGTGCTCTGCAAGAACAGTACAGTGTCTTAAAGATCCACTG
 GCTCATGGCAGCCTGGCGTTACGAAGAGTGTCTCCCTACTCTTCCACAGCATCACTACTACTTCATC
 AACAGCCAGGGCCACCCATCGAAGGCCCTCGCGTCATGCATTACATCACTCACCTGCTGAAGGGCGCCC
 TTCTCTTCATCACCATTGCCCTGATCGGCTCGGCTGGGCTTTGTGAAGTATATGCTGTGCGACAAGGA
 GAAGAAGATTTTTGGGATTGTGATTCCCTGCAGGTCTGGCTAATGTAGCATACATTGTATCGAGTCC
 CGTGAGGAGGGCGCCAGCGACTATGGAAGTCTGGAAGGAGATCTTGTTCCTGGTGGATCTCATCTGCTGTG
 GTGCCATCTCTCCAGTGGTCTGGTCCATCCGGCATCTACAGGACGCGTCTGGCACTGACGGGAAGGT
 GGCAATGAACCTGGCAAGCTGAAGCTGTTCCGACATTACTACGTATGGTCTGCTATATCTACTTT
 ACGCGCATCATGCCATTCTGCTCGAGTGGCGGTGCCCTTCCAGTGGCAGTGGCTGTACCAGCTCTTGG
 TGGAGAGTTCCACCCTGGCCTTCTCGTCTCACCAGTACAAGTCCAGCCGGCGGGGACAACCCATA
 CCTGCAGTGCACAGCAGGAGGATGAGGAGGACGTACAGATGGAGCAAGTATGACAGACTCTGGGTTT
 CGGAAGGCTGTCAAAGTTAACAAGACAGCCAGCGCCGGGAGCTGTTG

AC**CGGGCCGCT**CGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGA
 TTACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:

>RR207396 representing NM_199399
 Red=Cloning site Green=Tags(s)


MAVSERRGLSGESPAQCRWEYLSLLVLMVSGSGRIHRLTLTGEKRADIQLNSFGFYTNLSLEVELSLLR
 LSLQETEDKFPKVGFSLSRVRSVRSYSRNSHECPLERNSSNFLVFLINIKDLQVQVRKYGEQKLF
 SPGLLPEAPSSQSPPKPDPTGTPKDNHVIHPSPKMSAVKEDQAKLTVPQVSGDKALPAGHRHSSDQGPQ
 SQPPTRGPSGKEKDLVLGLGHLNDSYNFSFHIVIGSRAEEGQYSLNFHNCYNTIPQEQPFDLTMIREK
 NPEGLSAAEIPFLKLYLIMSACFLAAGIFWVSVLCKNTYSVFKIHWLMAALFTKSVSLLFHSINYYFI
 NSQGHPIEGLAVMHYITHLLKGALLFITIALIGSGWAFVKYMLSDKEKKIFGIVIPQLVLANVAYIVIES
 REEGASDYGLWKEILFLVDLICGAILFPVWVSIRHLQDASGTDGKVAMNLAKLFRHYVVMICYIYF
 TRIIAILLRVAVPFQWQWLYQLLVESSTLAFFVLTGYKFQPDGPNPYLQLPQQEDEDVQMEQVMTDSGF
 REGLSKVNKTASGRELL

TRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:

SgfI-NotI

Cloning Scheme:



ACCN: NM_199399

ORF Size: 1731 bp

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).


Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_199399.1](#), [NP_955431.1](#)

RefSeq Size: 2134 bp
 RefSeq ORF: 1734 bp
 Locus ID: 316136
 UniProt ID: [Q6P6V6](#)
 Cytogenetics: 9q12
 MW: 64.7 kDa

Product images:



Circular map for RR207396