

## Product datasheet for **RR207413L3V**

### **Pf4 (NM\_001007729) Rat Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	Pf4 (NM_001007729) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Pf4
Synonyms:	Cxcl4; PF-4; Pf4a; RATPF4A
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001007729
ORF Size:	315 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR207413).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001007729.1</a> , <a href="#">NP_001007730.1</a>
RefSeq Size:	318 bp
RefSeq ORF:	318 bp
Locus ID:	360918
UniProt ID:	<a href="#">P06765</a>
Cytogenetics:	14p22
Gene Summary:	member of low-molecular weight chemokines superfamily; may mediate inflammation and immune response [RGD, Feb 2006]



[View online »](#)