

## Product datasheet for **RR207453L4V**

### Hax1 (NM\_181627) Rat Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Hax1 (NM_181627) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Hax1
Synonyms:	HAX-1; Hs1bp1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_181627
ORF Size:	834 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR207453).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_181627.2</a> , <a href="#">NP_853658.1</a>
RefSeq Size:	1096 bp
RefSeq ORF:	837 bp
Locus ID:	291202
UniProt ID:	<a href="#">Q7TSE9</a>
Cytogenetics:	2q34
Gene Summary:	protein that may bind to bile salt export protein and may regulate its abundance in membranes [RGD, Feb 2006]



[View online »](#)