

Product datasheet for **RR207528**

Tspan1 (NM_001004236) Rat Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Tspan1 (NM_001004236) Rat Tagged ORF Clone
Tag: Myc-DDK
Symbol: Tspan1
Synonyms: MGC93753
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >RR207528 representing NM_001004236
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGCAGTGCTTCAAATTCATTAAGGTCATGATGATCCTCTTCAATCTACTCATCTTTCTCTGTGGTGCAG
 CCCTGTTGGCTGTGGGAATCTGGGTGTCCGTGACGGGACATCTTCTGAAGGCCTTCGGATCACTATC
 ATCGAGTGCCATGCAGTTTGTCAACGTGGGCTACTTCTCATCGCCGCTGGTGTACTCTTTATTCTT
 GGTTTCTGGGCTGCTATGGTGTCACTCCGAGAACAAGTGTGTGCTCATGATGTTCTTTCCATTCTCC
 TCATCATCTTCATTGCGGAGATTGCAGGGCCGTGGTTGCTTTGGTGTACACCACAATGGCTGAACAATT
 CCTGACATTTCTGGTGGTGCCTGCTATCGAAAAGGACTATGGTTACCAGACTGAATTCACCCAAGTGTGG
 AACTCTACGATGGAGGGGCTGCATTGCTGTGGCTTCAACAATTACACAGATTTAACTCCTCACGTTTCG
 TCAAAGAGAACAAGTCTTCTCCGTTCTGTTGCGCCAACAACACTGACAGTCATACAGTTGAACCATG
 CACCGAGGACAAGGCCAAAAGTATGAATGTACAGGGTTGTTTCAAACAGATTCTGCAGAAAATCAGAACC
 AATGCAGTTACCGTGGTGGTGTGGCAGTTGGAGTTGCGGCCCTAGAGCTGGCTGCCATGGTGGTATCCA
 TGTATCTATACTGCAATCTGAAG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >RR207528 representing NM_001004236
Red=Cloning site Green=Tags(s)

MQCFKFIKVMILFNLLIFLCGAALLAVGIWVSVDGTSFLKAFGSLSSSAMQFVNVGYFLIAAGAVLFIL
 GFLGCGYGAHSEKCVLMMFFSILLIIFIAEIAAGAVVALVYTTMAEQFLTFLVVPATIEKDYGYQTEFTQVW
 NSTMEGLHCCGFNNTDFNSSRFVKENKVFPPFCCANNTDSHTVEPCTEDKAKSMNVQGCQKILQKIRT
 NAVTVGGVAVGVAELEAAMVSMYLYCNLK

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/ja2447_f01.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF

ACCN: NM_001004236

ORF Size: 723 bp

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001004236.1](#), [NP_001004236.1](#)

RefSeq Size: 1256 bp

RefSeq ORF: 726 bp

Locus ID: 298436


UniProt ID: [Q6AYR9](#)

Cytogenetics: 5q35

MW: 26.9 kDa

Gene Summary: human homolog is associated with cervical neoplastic cell proliferation and cervical carcinoma [RGD, Feb 2006]

Product images:



Circular map for RR207528