

Product datasheet for **RR207562**

Aspa (NM_024399) Rat Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Aspa (NM_024399) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Aspa
Synonyms:	MGC93407; MGC93408
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin
ORF Nucleotide Sequence:	>RR207562 representing NM_024399 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGACTTCTTGTGTGGCTGAAGAACCTATTA AAAAGATTGCCATCTTCGGAGGGACTCATGAAATGAAC
TGACTGGAGTGTCTAGTTACTCACTGGCTAAAGAATGGCGCTGAAGTTCACAGAGCAGGGCTGGAAGT
GAAGCCATTCATTACCAACCAAGAGCGGTGGAGAAGTGCACCAGATACATTGACTGTGACCTGAACCGT
GTTTTTGACCTTGAAAATCTTAGCAAAGAGATGTCTGAAGATTTGCCGTATGAAGTTAGAAGGGCTCAAG
AAATAATCACTATTTGGTCCAAAAGATAGTGATGATGCCTACGATGTCGTTTTTGTCCCTCATAACAC
TACTTCTAACATGGGTTGCACTCTTATTCTTGAGGATTCCAGGAATGACTTTTTAATCCAGATGTTTCAC
TATATTAAGACCTGCATGGCTCCATTACCTGCTCTGTTACCTCATCGAGCATCCTTCCCTCAAATATG
CAACCACTCGTTCCATTGCCAAGTATCCTGTTGGTATAGAAGTTGGTCTCAGCATCACGGTGTCTGAG
AGCTGATATTTTAGACCAAATGAGAAGAATGCTCAAGCAGCTCTTGATTTTATACAGCGTTTCAATGAA
GGTAAAGAATTTCCCCCTGTGCTATTGATGTCTATAAAAATAATGGAGAAAAGTTGATTATCCAAGGAACG
AAAGTGGAGACGTGGCTGCTGTTATCCATCCTAATCTGCAGTGCAGGATCAAGACTGGAACCAATTGCA
CCCTGGAGATCCTGTGTTGTGCTCTTGATGGAAAAGTTATTCCACTGGGTGGAGACTGTACCGTGTAC
CCAGTGTGTTGTAATGAAGCTGCATATTATGAAAAAAGAAGCATTTGCAAAGACAACAAAACCTAACAC
TCAACGCAAAAAGCATCCGCTCCACTTTGCAC

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_024399.1 , NP_077375.1
RefSeq Size:	1552 bp
RefSeq ORF:	945 bp
Locus ID:	79251
UniProt ID:	Q9R1T5
Cytogenetics:	10q24
MW:	35.6 kDa
Gene Summary:	catalyzes the deacetylation of N-acetylaspartate to free acetate; may play a role in central nervous system myelination [RGD, Feb 2006]