

Product datasheet for **RR207593**

Trappc6b (NM_001106733) Rat Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Trappc6b (NM_001106733) Rat Tagged ORF Clone
Tag: Myc-DDK
Symbol: Trappc6b
Synonyms: RGD1309325
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RR207593 representing NM_001106733
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGGCGGACGAGGCGTTGTTTTGCTTCTCCATAACGAGATGGTGTCCGGAGTATATAAGTCCGCCGAGC
AGGGGGAAGTGGAAAATGGACGGTGTACTAAGCTGGAAAACATGGGGTTTCGAGTGGGACAAGGATT
GATAGAAAAGTTTACAAAAGATACTGCAAGTTCAAGGATGAATTAGACATCATGAAGTTCATATGTAAA
GATTTTTGGACTACGGTATTCAAGAAACAAATTGACAATCTAAGGACAAATCATCAGGGCATATATGTGC
TTCAGGACAACAAATTCGACTACTCATTCAAGTGTCTGCAGGAAAACAGTATTTAGAACATGCATCCAA
GTATCTAGCGTTTACATGTGGCTTAATCAGAGGTGGCTTGTCAAACCTGGGAATAAAAAGTATTGTGACA
GCTGAAGTATCTTCAATGCCTGCCTGCAAATCCAGGTGATGATACAGAAGCTG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RR207593 representing NM_001106733
Red=Cloning site Green=Tags(s)


MADEALFLLLNEMVSGVYKSAEQGEVENGRICITKLENMGFRVGGQGLIERFTKDTARFKDELDMKFICK
DFWTTVFKKQIDNLRNTHQGIYVLQDNKFRLLIQLSAGKQYLEHASKYLAFTCGLIRGGLSNLGIKSIIVT
AEVSSMPACKFQVMIQKL

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI



Cloning Scheme:

Plasmid Map:


ACCN: NM_001106733

ORF Size: 474 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001106733.1, NP_001100203.1</u>
RefSeq Size:	1165 bp
RefSeq ORF:	477 bp
Locus ID:	299075
Cytogenetics:	6q23
MW:	18 kDa