

## Product datasheet for **RR207602L4V**

### Cwc27 (NM\_001013199) Rat Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Cwc27 (NM_001013199) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Cwc27
Synonyms:	Sdccag10
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001013199
ORF Size:	1404 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR207602).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001013199.1</a> , <a href="#">NP_001013217.1</a>
RefSeq Size:	1852 bp
RefSeq ORF:	1407 bp
Locus ID:	361887
UniProt ID:	<a href="#">Q5XIB2</a>
Cytogenetics:	2q12
Gene Summary:	As part of the spliceosome, plays a role in pre-mRNA splicing. Probable inactive PPlase with no peptidyl-prolyl cis-trans isomerase activity.[UniProtKB/Swiss-Prot Function]



[View online »](#)