

Product datasheet for **RR207609**

Glt8d1 (NM_001007683) Rat Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Glt8d1 (NM_001007683) Rat Tagged ORF Clone
Tag: Myc-DDK
Symbol: Glt8d1
Synonyms: Da2-24
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RR207609 representing NM_001007683
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGTCATTCCGTAAGTAACCATCATCATCTGGGCCCTGGCCGTCATTCTTCTACTGGCTTTGCACC
ATAACTTCCTCAGCTTGAGCAGTTTGTTAAGGAATGACGTTTCAGATTCAGGAATTGTGGGGCTTCAGCC
CATAGACTTTACCGCAAATTTTCATCAACATCCAGTAAATGAGAGACAAGAAGAGATTCGGTGGTTATT
GCTGCCTCTGAAGATAGGCTGGGTGGGACCATTGCAGCCATAAATAGCATTTCATAAAACTCGCTCTA
ATGTAATTTTTTACATTGTTACATTCAATAGAACAGCGGACCATCTCAGGTCCTGGCTTAACAGTGGTTC
CCTGAAAAGCATCAGGTACAAAATTGTAATTTTGACACTAAACTTTTGGAAAGGAAAGTAAAGGAGGAT
CCTGACCAGGGGGAGTCCATGAAACCATTAACCTTTGCAAGGTTCTACTTGCCAATTCTGGTTCAGTG
CAAAGAAGGCCATTTACATGGATGACGATGTAATTGTACAAGGTGATATACTTGCCTTATAACAATACCC
ACTGAAGCCAGGACATGCTGCTGCATTTTCAGAAGACTGTGACTCAGCTTCTACTAAAGTCATGATCCGT
GGAGCAGGGAACCAAGTACAATTACATTGGTTATCTTGACTATAAAAAGGAAAGAATTCGTAAGCTATCCA
TGAAAGCCAGCACCTGCTCCTTTAATCCTGGAGTTTTTGTGCAACTTGACTGAGTGGAAAAGACAGAA
TGTAACTAACCAGCTGGAAAAATGGATGAAGCTCAATGTAGAAGAGGGACTGTACAGCAGAACCCTGGCT
GGCAGCATCACACCTCCTCTGCTTATCGTATTTTATCAGCAGCATTCCACCATTGACCTATGTGGA
ATGTCCGCCACCTTGGTTCAGTCTGGAAAACGGTATTCACCCAGTTTGTAAGGCTGCCAAGTTACT
TCACTGGAATGGTCACTTCAAGCCTTGGGAAGAGCTGCTTCATATGCTGACGTTTGGGAGAAATGGTAT
ATTCCAGACCCAACAGGCAAATTTAGTTAATCCGAAGACATATGGATACCTCAAACATAAAG

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >RR207609 representing NM_001007683
Red=Cloning site Green=Tags(s)


MSFRKVTIIWALAVILFLLALHHNFLSLSSLLRNDVSDSGIVGLQPIDFTANFHQHPVNERQEEIPVVI
 AASEDRLLGGTIAAINS IHQNTRS NVIFYIVTFNRTADHLRSWLN SGLSKSIRYKIVNFDTKLLEGKVKED
 PDQGESMKPLTFARFYLPILVPSAKKAIYMDDDVIVQGDILALYNTPLKPGHAAAFSEDCDSASTKVMIR
 GAGNQYNYIGYLDYKKERIRKLSMKASTCSFNPGVFVANL TEWKRQNVTNQLEKWMKLNVEEGLYSRTLA
 GSITTPLLIVFYQQHSTIDPMWNVRHLGSSAGKRYSPQFVKAAKLLHWNHGFKPWGAAASYADVWEKWY
 IPDPTGKFSLIRRHMDTSNIK

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Restriction Sites:

Sgfl-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001007683

ORF Size:	1113 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001007683.1 , NP_001007684.1
RefSeq Size:	1644 bp
RefSeq ORF:	1116 bp
Locus ID:	306253
UniProt ID:	Q6AYF6
Cytogenetics:	16p16
MW:	42.1 kDa