

## Product datasheet for RR207634L4

### Inpp4b (NM\_053917) Rat Tagged Lenti ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Inpp4b (NM_053917) Rat Tagged Lenti ORF Clone
Tag:	mGFP
Symbol:	Inpp4b
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR207634).
Restriction Sites:	Sgfl-MluI
Cloning Scheme:	

Cloning sites used for ORF Shuttling:



\* The last codon before the Stop codon of the ORF.



[View online »](#)



RefSeq ORF: 2829 bp  
Locus ID: 116699  
UniProt ID: [Q9QWG5](#)  
Cytogenetics: 19q11  
Gene Summary: enzyme that catalyzes the hydrolysis of the 4-position phosphate of phosphatidylinositol 3,4-bisphosphate, inositol 1,3,4-trisphosphate, and inositol 3,4-bisphosphate [RGD, Feb 2006]