

Product datasheet for **RR207637L3V**

Asic1 (NM_024154) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Asic1 (NM_024154) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Asic1
Synonyms:	Accn2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_024154
ORF Size:	1578 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR207637).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_024154.2 , NP_077068.1
RefSeq Size:	3771 bp
RefSeq ORF:	1581 bp
Locus ID:	79123
UniProt ID:	P55926
Cytogenetics:	7q36
Gene Summary:	forms a proton-gated ion channel that may play a role in sensory neuron function and pain perception [RGD, Feb 2006]



[View online »](#)