

Product datasheet for **RR207659**

Fbxo28 (NM_001107203) Rat Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Fbxo28 (NM_001107203) Rat Tagged ORF Clone
Tag: Myc-DDK
Symbol: Fbxo28
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RR207659 representing NM_001107203
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGCGGGCGGCTCTGAGGAGCGGATGGCTGAGGAAGGAGGGCGGCCACGGCGACGGCGGCTCCTGT
CGGCCCGCAGCTCCGCTCAGCGACAACCCCGACACCCCGTCGACGGCCCTGCAGCCGGGTCTCAGGC
ACCCGCCGCGCCGCGCTGGCACCGGACCACCTGCCTCAAACAACACACTGGTGGCGCTGCCCATCGTA
GCCATCGAGAACATCCTCAGCTTTATGTCTACGACGAAATTAGCCAGCTCCGCCTGGTTTGTAAGAA
TGGACTTGGTCTGTGAGAGAATGCTGAATCAGGGCTTTCTGAAAGTGGAGAGGTACCATAACCTGTGTC
GAAACAAGTTAAAGCACAACTACCAAGGAGGGAGTCAGAAAGGAGGAACCACTCTTAGCTCGCCATGCA
GACATCCTTGCTGCTGTGAAACGCGGCTGTCACTGTTAAACATGACCTTCATGAAGTATGTGGACTCCA
ACCTCTGTTGCTCATCCAGGAAAGGTGATTGATGAGATTTACCGTGTACTGAGATATGTCAATTCTAC
CAGAGCCCTCAGCGAGCTCATGAAGTCTGCAGGAGTTGCGGGACATCTCCTCCATGGCAATGGAGTAC
TTTGATGAGAAGATTGTTCCGATCTTGAAGAGGAAATTGCCAGGATCCGATGTGTCTGGAAGACTCATGG
GCTCTCCTCCAGTCCCAGGACATCCGACGCCCTGACAACCATGCAGCTTCTCCAAGCAAACCCGTC
CAGACAAGAGGTTACAAACTGCAGCAGAGGTCAAAACCAATGGTGTGGCGTACTGTCTCAGGCGT
GAGATTTCTGAGCTCCGACCAAAGTGCAAGAACAGCAGAAACAGCTGCAAGACCAAGACCAAGAACTGC
TAGAGCAGACCCAGATCATAGGCGAGCAAAACGCTCGGCTGGCGGAGCTGGAGCGGAAGCTCCGAGAAGT
CATGGAAGTGCTGTAGGAACATCCTCAGGGTCCGGGCAGAGCGAGGAGTCTCCTCGAAACGAAGGAAG
GCAACAGAAGCCATAGACTCACTTAGGAAATCAAACGACTGCGGAATAGGAAG

ACGCGTACGCGGGCGGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)


Protein Sequence: >RR207659 representing NM_001107203
 Red=Cloning site Green=Tags(s)

MAAASEERMAEEGGGGHGDGGSCSAASSAQRQPPTPPSQALQPGSQAPAAPALAPDHL PQNNTLVALPIV
 AIENILSFMSYDEISQLRLVCKRMDLVCQRMLNQGFLKVERYHNL CQKQVKAQLPRRESERRNHSLARHA
 DILA AVETRL SLLNMTFMKYVDSNLCCFIPGKVIDE IYRVLRYVNSTRAPQRAHEVLQELRDISSMAMEY
 FDEKIVPILKRRKLPGSDVSGRLMGSPVPGPSAALTTMQLFSKQNP SRQEVTKLQQQVKTNGAGVTVLRR
 EISELR TKVQEQQKQLQDQDKLLEQTQI IGEQNARLAELERKLRVME SAVGTSSGSGQSEESPRKRRK
 ATEAIDSLRKSRLRN RK


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001107203

ORF Size:	1104 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001107203.3, NP_001100673.1</u>
RefSeq Size:	1523 bp
RefSeq ORF:	1107 bp
Locus ID:	305105
Cytogenetics:	13q26
MW:	41.1 kDa