

Product datasheet for **RR207684L4V**

Pdlim3 (NM_053650) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Pdlim3 (NM_053650) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Pdlim3
Synonyms:	Actn2lp; Alp
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_053650
ORF Size:	1086 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR207684).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_053650.1 , NP_446102.1
RefSeq Size:	1568 bp
RefSeq ORF:	1089 bp
Locus ID:	114108
UniProt ID:	Q66HS7
Cytogenetics:	16q11
Gene Summary:	binds to the spectrin-like motifs of alpha-actinin-2 in skeletal muscle [RGD, Feb 2006]



[View online »](#)