

Product datasheet for **RR207695**

Cyfp2 (NM_001106996) Rat Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Cyfp2 (NM_001106996) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Cyfp2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RR207695 representing NM_001106996
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGC**

ATGTCCTCTACCCCTTGATCTGTACAACGACAGCGCCTACTATGCGCTACCAAGTTCAAAAAGCAGTT
 CCTGTATGACGAGATCGAAGCAGAGGCAAGTGAACCTGTGCTTTGATCAGTTTGTCTACAAGCTGGCTGA
 CCAGATCTTTGCCTACTACAAAGCCATGGCCGGCAGTGTCTGTTGGACAAACGTTTTAGAGCTGAATGT
 AAGAACTATGGAGTGATCATTCCGTACCGCCATCTAACCGGTACGAGACACTGCTGAAACAGAGACACG
 TCCAGCTGTTGGGTAGATCAATTGACTTGAACAGACTCATCACGACGCGGATCTCTGCCGCCATGTACAA
 ATCCTTGGACCAGGCAATCAGCCGCTTTGAGAGCGAAGACCTGACCTCCATTGTGGAGCTAGAGTGGCTT
 CTGGAGATTAACCGGCTCACACACCGACTCCTCTGTAAGCACATGACGCTAGACAGCTTCGACGCCATGT
 TCCGAGAAGCCAATCATAATGTGTCTGCCCCCTATGGCCGCATCACCTGCATGTCTTCTGGGAGCTGAA
 CTTTGACTTCTCCCAACTACTGCTATAACGGCTCCACCAACCGCTTTGTCCGCACTGCCATCCCTTTC
 ACCCAAGAACCACAACGGGATAAGCCCGCCAATGTCCAGCCTTACTACCTCTACGGATCCAAGCCTCTCA
 ACATTGCCTACAGCCACATCTACAGCTCCTACCGCAACTTTGTGGGGCCCCCTCACTTCAAGACCATCTG
 CAGACTGCTGGGTACCAGGGCATCGCCGTGGTCATGGAGGAAGTCTGAAAAATCGTCAAGAGCTTGCTC
 CAGGGTACCATTCTGCAGTACGTGAAGACACTCATAGAGGTGATGCCTAAAAATGCCGTTTGCCTCGTC
 ACGAGTATGGCTCTCCAGGGATCCTGGAGTTCTTCCACCACCAACTGAAGGACATCATCGAGTATGCAGA
 ACTCAAAACAGACGTGTTCCAGAGTCTCAGAGAGGTGGCAATGCCATCCTGTTCTGCCTCCTCATAGAG
 CAAGCTCTGTCTCAAGAGGAAGTCTGTGATTTGCTCCATGCCGCACCCTTCCAAAACATTCTGCCTCGGG
 TCTATATCAAAGAGGGGGAGCGTCTGGAAGTCCGGATGAAACGCCCTCGAGGCCAAGTACGCCCCACTTCA
 CCTGGTCCCCCTGATAGAGCGACTGGGGACGCCCTCAGCAAATTGCCATCGCTCGGGAGGGTGACCTGCTG
 ACCAAGGAGCGGCTCTGCTGTGGCTGTCCATGTTGAGGTTATCCTGACCCGAATTTCGGAGCTACCTGC
 AGGACCCCATCTGGAGGGGCCGCCACCCACCAATGGCGTCATGCACGTGGATGAGTGTGTTGAGTTTCA
 CCGGCTGTGGAGCGCCATGCAGTTTGTCTACTGCATCCCTGTGGGGACCAACGAGTTCACAGCTGAGCAG
 TGCTTCGGTGATGGCTTGAAGTGGGCTGGCTGCTCTATCGTTGTCTTGTGGGCCAGCAACGGCGCTTTG
 ACCTGTTTGATTTCTGTTATCATCTGCTGAAGGTGCAGAGGCAGGACGGGAAGGATGAGATCATTAAAGAA
 TGTGCCCTGAAGAAGATGGCTGACCGCATCAGGAAGTACCAGATCTTGAACAATGAGGTTTTTGCCATT
 CTGAACAAGTACATGAAGTCCGTAGAGACAGACAGCTCCACGGTAGAGCACGTACGCTGCTTCAGCCCC
 CAATCCACCAGTCCCTGGCCACCACCTGC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:

>RR207695 representing NM_001106996
 Red=Cloning site Green=Tags(s)


MSSTPWICTTTAPTMRSPSSKSSSCMTRSKQRQVNLCFDQFVYKLADQIFAYYKAMAGSVLLDKRFRAEC
 KNYGVIIIPYPPSNRYETLLKQRHVQLLGRSIDLNRLITQRISAAMYKSLDQAISRFESEDLTSIVELEWL
 LEINRLTHRLLCKHMTLDSFDAMFREANHNVSAPYGRITLHVFWELNFDLPNYCYNGSTNRFVRTAIPF
 TQEPQRDKPANVQPYLYGSKPLNIAYSHIYSSYRNFGVPPHFKTICRLLGYQGIADVMEELLKIVKSL
 QGTILQYVKTLIEVMPKICRLPRHEYSGPILFEFFHQLKDIIEYAEKTDVVFQSLREVGNAILFCLLIE
 QALSQEEVCDLLHAAPFQNILPRVYIKERLEVRMKRLEAKYAPLHLVPLIERLGTTPQQIAIAREGDL
 TKERLCCGLSMFEVILTRIRSYLDPIWRGPPPTNGVMHVDECVEFHRLWSAMQFVYCIPVGTNEFTAEQ
 CFGDGLNWAGCSIVVLLGQRRFDLDFCYHLLKVQRQDGKDEIKNVPLKKMADRIRKYQILNNEVFAI
 LNKYMKSVETDSSSTVEHVRCFQPPIHQSLATTC

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:

Sgfl-MluI

Cloning Scheme:



ACCN: NM_001106996

ORF Size: 1779 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001106996.1](#), [NP_001100466.1](#)

RefSeq Size: 6253 bp


RefSeq ORF: 1782 bp

Locus ID: 303073

Cytogenetics: 10q21

MW: 68.7 kDa

Product images:



Circular map for RR207695