

Product datasheet for RR207719L4

Pou2f3 (NM_001105745) Rat Tagged Lenti ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Pou2f3 (NM_001105745) Rat Tagged Lenti ORF Clone
Tag:	mGFP
Symbol:	Pou2f3
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR207719).
Restriction Sites:	Sgfl-MluI
Cloning Scheme:	

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF.



[View online »](#)

RefSeq ORF: 1293 bp
Locus ID: 116544
Cytogenetics: 8q22
Gene Summary: may play a role in epidermal development; variant Skn-1i inhibits Oct1 transactivation; variant Skn1a activates cytokeratin 10 (K10) gene expression [RGD, Feb 2006]