

Product datasheet for **RR207725L4V**

Hdgf (NM_053707) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Hdgf (NM_053707) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Hdgf
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_053707
ORF Size:	711 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR207725).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_053707.2 , NP_446159.1
RefSeq Size:	1847 bp
RefSeq ORF:	714 bp
Locus ID:	114499
UniProt ID:	Q8VHK7
Cytogenetics:	2q34
Gene Summary:	induces DNA replication in vascular smooth muscle cells [RGD, Feb 2006]



[View online »](#)