

Product datasheet for **RR207734**

Snrpa (NM_001008303) Rat Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Snrpa (NM_001008303) Rat Tagged ORF Clone
Tag: Myc-DDK
Symbol: Snrpa
Synonyms: MGC105915
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RR207734 representing NM_001008303
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGGCAGTCCCGAGACCCGTCCCAACCACACTATTTATCAACAATCTCAATGAGAAGATCAAGAAAG
ATGAGCTCAAGAAGTCCCTGTATGCCATCTTCTCCAGTTTGCCAGATCCTGGATATCCTGGTGTCTCG
GAGCCTGAAGATGAGGGCCAGGCCTTCGTTCATCTTCAAGGAGGTCAGCAGTGCCACCAACGCCCTGCGT
TCCATGCAGGGCTTCCCTTTCTACGACAAGCCATGCGTATCCAGTACGCAAGACTGACTCAGACATCA
TTGCCAAAATGAAGGGCACCTTTGTGGAGAGAGATCGAAACGTGAGAAGAGGAAGCCCAAGAGTCAGGA
GACACCTGCTGCCAAAAGGCTGTTTCAGGGTGGGGCAGCTGCCCTGTGGTTGGCGCTGTCCAGCCTGTA
CCGGGCATGCCACCGATGACTCAGGCACCCCGCATCATGCACCACATGCCAGGACAGCCTCCCTACATGC
CACCACCTGGCATGATCCCGCCACCAGGCCTTGTCTCGCCAGATCCCTCCTGGGGCAATGCCCCACA
GCAGCTCATGCCTGGGCAGATGCCACCTGCCAACCTCTCTCGGAGAATCCACCAATCACATCCTGTTT
CTCACCAACCTGCCTGAGGAGACCAACGAGCTCATGCTCTCCATGCTCTCAACCAGTCCCTGGCTTCA
AGGAGGTGCGTCTGGTCCCTGGGCGCCATGACATTGCCTTCGTGGAGTTTGACAATGAAGTGCAGGCTGG
GGCAGCAGAGATGCCCTGCAAGGCTTCAAGATCACACAGAACAATGCAATGAAGATCTCTTTTGCCAAG
AAG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >RR207734 representing NM_001008303
Red=Cloning site Green=Tags(s)


MAVPETRPNHTIYINNLNEKIKKDELKKSLEYAIFSQFGQILDILVSRSLKMRGQAFVIFKEVSSATNALR
 SMQGFPFYDKPMRIQYAKTDSIDI AKMKGTFFVERDRKREKRKPKSQETPAAKKAVQGAAAPVVGAVQPV
 PGMPPMTQAPRIMHHMPGQPPYMPPIPPGLAPGQIPPGAMPQQQLMPGQMPPAQPLSENPPNHILF
 LTNLPEETNELMLSMLFNQFPGFKEVRLVPRHDI AFVEFDNEVQAGAARDALQGFKITQNNAMKISFAK
 K

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Restriction Sites:

Sgfl-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001008303

ORF Size: 843 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001008303.1 , NP_001008304.1
RefSeq Size:	1170 bp
RefSeq ORF:	846 bp
Locus ID:	292729
Cytogenetics:	1q21
MW:	31.2 kDa