

Product datasheet for **RR207757**

Taf15 (NM_001105824) Rat Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Taf15 (NM_001105824) Rat Tagged ORF Clone
Tag: Myc-DDK
Symbol: Taf15
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >RR207757 representing NM_001105824
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGAAATCAAGGCAGTCAAGGCTATGGACAGACACAACAAAGCTATTCTGGCTATGGGCAAACGACTGA
TTCTTCATATGGACAAAATATGGTGGTTACTCCGGTTATGGACAAAATCAATCAGGTTATTCACAGTCC
TATGGTAGTTATGAGAATCAAAGCAGAGTTCCTATGGCCAGCAGTCATATAATAACCAGGGACAACAAA
ACACGGAATCATCAGGAGGTCAAGGTGGAAGAGCACCTTCATATGGCCAGTCAGACTATGGTCAGCAAGA
TTCTTATGACCAGCAGTCAGGTTATGATCAGCACCAAGGCTCATATGATGAGCAGTCCAATTATCAGCAG
CATGATTCCTATAATCAAACCAGCAGTCTTACCATCCACAAAGGAAAACTATAGCCACCACACACAAG
TGCAGTGGTGGTGACCGGGTGGCTTCAAAAATTTGGTGGTCACAGGGATTATGGACCCAGACCAGATG
CTGATTCAGAGTCTGATAATTCAGATAACAACACAATCTTTGTGCAAGGACTTGGGGAGGGTGTGTCTAC
AGATCAAGTGGGAGAGTTCTTTAAGCAAAATAGGAATTATCAAGACTAACAAAAAGACTGGAAAACCAATG
ATAAATCTTTACACAGACAAAGACACAGGAAAGCCAAAAGGAGAAGCAACAGTGTCAATTTGATGACCCTC
CATCAGCTAAGGCAGCCATTGATTGGTTTGTGAAAAGAATTCATGAAAACATCATTAAAGTGTCAAT
TGCCACGAGAAGGCCTGAATTCATGAGAGCGGTGGAAGTGGAGGTGGACGGCGAGGTCTGGAGGATAT
AGAGGCCGTGGAGTTTTAGGGGCGAGGTGGAGACCCCAAAAATGGGGACTGGGTTTGCCCAATCCGT
CATGTGGGAATATGAATTTGCTCGAAGGAATTCCTGTAATCAGTGCAATGAGCCAAGACCAGGAGACTC
TCGTCCCTCAGGGGAGAAATGACTACAGAAATGATCAGCGCAACCGACCGTACTGATGACTGTTTCCAA
TACTCCCTGTCTCTGACATGATCCATGGGAATTGTCAGAGTGTGCCTGCTGCTTTCCTCATGGCCTCT
TGAGTAGTGAAATTGACATTTGGTTTTATTTGGGTGGGAGGGCTGGGACAGTTTTCTT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)


Protein Sequence: >RR207757 representing NM_001105824
Red=Cloning site Green=Tags(s)

```
MEIKAVKAMDRHNKAILAMGKRLILHMDKTMVVTPTVMDKINQVIHSPMVMRIKSRVPMASSHIITRDNK
TRNHQEVKVEEHLHMASQTMVSKILMTSSQVMISTKAHMMSSPIISSMIPIIKTSSLTIHKGKTIATTHK
CSGGDRGGFKNFGGHRDYGPRPDADSESDNSDNTIFVQGLGEGVSTDQVGEFFKQIGIIKTKKTKGKPM
INLYTDKDTGKPKGEATVSFDDPPSAKAAIDWFDGKEFHGNIKIVSFATRRPEFMRGGSGGRRRGRGGY
RGRGGFQGRGGDPKNGDWVCPNPSCGNMNFARRNSCNQCNEPRPESRPSGGEMTTEMISATDRTDCCFQ
YSLVSDMIHGELSECACCFPHGLLSSEIDIWVFIWVGGLGTVFL
```

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001105824

ORF Size: 1182 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001105824.1](#), [NP_001099294.1](#)

RefSeq Size: 1494 bp


RefSeq ORF: 1185 bp

Locus ID: 287571

Cytogenetics: 10q26

MW: 43.2 kDa

Product images:



Circular map for RR207757