

Product datasheet for **RR207775**

Unc5b (NM_022207) Rat Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Unc5b (NM_022207) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Unc5b
Synonyms:	Unc5h2
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

ORF Nucleotide
Sequence:

>RR207775 representing NM_022207
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCCGCATCGCC

ATGAGGGCCCGAGCGGGCGCGGGCGCGCTGCTGCTGGCGCTGCTGCTGCTGCTGGGATCCGACACCGA
 GCTTAGCAGGCATTGACTCTGGTGGCCAGGCACTCCAGACTCCTTCCATCAGCACCCGCGGAGCAGCT
 GCCTCACTTCTGCTGGAACAGAGGATGCCTACATCGTAAAGAACAAGCCAGTGAATGCACTGCCGA
 GCCTTCCCTGCCACACAGATCTACTTCAAGTGAATGGCGAGTGGGTTAGCCAGAAAGGCCACGTACGC
 AGGAGAGCCTGGATGAGGCCACAGGCTTGCGAATACGAGAGGTGCAGATAGAGGTGTGCGGCAGCAGGT
 GGAGGAACTTTTGGGCTCGAGGACTACTGGTGTGAGTGCCTGGGCTCTTCGGAAACCACCAAG
 AGTCGCCGAGCCTACATCCGATTGCCTACTTGCAGCAAGAACTTTGACCAGGAGCCTCTGGCGAAGGAGG
 TACCCTTGGATCATGAGGTCTTCTGCAGTCCCGCCACCAGAGGGAGTGCCTGTGGCTGAGGTGGAATG
 GCTCAAGAATGAAGATGTCATCGATCCCGCTCAGGACACTAACTTCTGCTCACCATTGACCACAACCTC
 ATCATCCGCCAGGCGCGCTCTCAGACACAGCCAACTACACCTGTGTGGCAAAGAATATTGTGGCCAAGC
 GCCGGAGCACGACGGCCACAGTCATCGTCTATGTGAACGGAGGTTGGTCCAGCTGGGCAGAAATGGTACC
 CTGCTCTAACCGCTGCGGCCGAGGTTGGCAGAAACGTAAGTACTAGGACCTGCACCAACCCAGCCCACTCAAT
 GGAGGTGCCTTCTGCGAGGGACAGGCTTGCAGAAAGACGGCTTGCACCACCGTGTGCCAGTGGATGGAG
 CGTGGACTGAGTGGAGCAAGTGGTCCGCTGCAGCACAGAGTGTGCGCACTGGCGCAGCCGCGAGTGCAT
 GGCACCCGCCCCAGAACCGAGGCCGTGACTGCAGCGGGACGCTACTTACTCCAAGAAGTGCACCGAT
 GGGCTGTGCGTGTGAATCAGAGAACTTAAACGACCTAAAAGCCGCCCCCTGGAGCCGTCGGGAGACG
 TGGCGTGTATGCGGGCCTCGTGGTGGCGTCTTGTGGTCTGCGAGTCTCATGGCTAGGAGTGTGAT
 CGTGTACCGGAGAACTGCCGGACTTCGACACGGACATCACTGACTCCTCTGCTGCCCTCACTGGTGGT
 TTCACCCCGTCAACTTCAAGACTGCAAGGCCAGCAACCCACAGCTCCTGCACCCATCCGCCCTCCGG
 ACCTAACGGCCAGTGTGGCATCTACCGCGACCTGTGTATGCCCTGCAGGACTCTGCCACAAGATCCC
 TATGACTAATCACCCCTTCTGGATCCCTTGCCAGCCTCAAGATCAAGGTCTATGACTCCAGCACCATC
 GGCTCTGGGGCTGGCCTGGCTGATGGAGCCGACCTGCTGGGTGTCTTACCACCCGGTACATACCCAGGCG
 ATTTCTCCGGGACACCCACTTCTGCACCTGCGCAGCGCCAGCCTTGGTCCCAGCACCTCTGGGCT
 CCCTCGAGACCCAGCAGCAGTGTGAGTGGCACCTTGGTGGCTGGGTGGGAGGCTGACCATCCCGGC
 ACAGGGGTGAGCCTGTTGTTACCAATGGAGCCATTTCCAGGGCAAGTCTATGACTTGTATCTACGTA
 TCAACAAGACTGAAAGCACCCCTCCCACTTTCGGAAGGTTCCAGAGCAGTATTGAGCCCTCGGTGACCTG
 CGGGCCACGGGCTCCTCCTGTGCCGCTGTTGTCCTCACTGTGCCCACTGTGCTGAAGTATTGCC
 GGAGACTGGATCTTCCAGCTCAAGACCCAGGCCATCAGGGCCACTGGGAGGAGGTGGTACTTTGGATG
 AGGAGACTCTGAACACCCCTGCTACTGCCAGCTAGAGGCTAAATCCTGCCACATCCTGTTGGACCAGCT
 GGGTACCTACGTGTTACGGGCGAGTCTACTCCCGCTCCGCAAGTCAAGCGGCTCCAGCTAGCCATCTTC
 GCCCAGCCCTGACCTCCCTGGAGTATAGTCTCAGGGTCTACTGTCTGGAGGACACTCTGCAGCAC
 TGAAGGAGGTCTAGAGCTGGAGAGGACTCTGGGTGGCTACTTGGTGGAGGAGCCCAAGACTTTGCTCTT
 TAAGGACAGTTACCACAACCTACGNCTCCCTCCATGACATCCCCATGCCCACTGGAGGAGCAAACTA
 CTGGCCAAGTACCAGGAGATTCCTTCTACCATGTGTGGAACGGCAGCCAGAAAGCCCTGCACTGCACTT
 TCAACCCTGGAGAGACATAGCCTAGCCTCACTGAGTTCACCTGTAAGGTCTGCGTGCAGGAGTGAAGG
 GGAAGGCCAGATTTTCCAGCTGCACACCAGCTGGCTGAGACGCCTGCTGGCTCCCTGGATGCACTCTGC
 TCTGCCCTGGCAATGCTGCCACCACACAGCTGGGACCCATGCTTCAAGATACCAGTGTCCATCCGCC
 AGAAGATCTGCAACAGCCTGGACGCCCAACTCACGGGCAATGACTGGCGGCTGTTGGCACAGAAGCT
 CTCCATGGACCGTACCTGAAGTACTTCGCCACCAAGCTAGTCCACAGGCGTGTCTTAGACCTCTGG
 GAAGCTCGGCAGCAGGATGATGGGACCTCAACAGCCTGGCCAGTGCCTTGGAGGAGATGGCAAGAGTG
 AGATGCTGGTAGCCATGACCACTGATGGCGATTGC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RR207775 representing NM_022207
 Red=Cloning site Green=Tags(s)

MRARSGARGALLALLLWDP T P S L A G I D S G G Q A L P D S F P S A P A E Q L P H F L L E P E D A Y I V K N K P V E L H C R
 A F P A T Q I Y F K C N G E W S Q K G H V T Q E S L D E A T G L R I R E V Q I E V S R Q Q V E E L F G L E D Y W C Q C V A W S S S G T T K
 S R R A Y I R I A Y L R K N F D Q E P L A K E V P L D H E V L L Q C R P P E G V P V A E V E W L K N E D V I D P A Q D T N F L L T I D H N L
 I I R Q A R L S D T A N Y T C V A K N I V A K R R S T T A T V I V Y V N G G W S S W A E W S P C S N R C G R G W Q K R T R T C T N P A P L N
 G G A F C E G Q A C Q K T A C T T V C P V D G A W T E W S K W S A C S T E C A H W R S R E C M A P P P Q N G G R D C S G T L L D S K N C T D
 G L C V L N Q R T L N D P K S R P L E P S G D V A L Y A G L V V A V F V L A V L M A V G V I V Y R R N C R D F D T D I T D S S A A L T G G
 F H P V N F K T A R P S N P Q L L H P S A P P D L T A S A G I Y R G P V Y A L Q D S A D K I P M T N S P L L D P L P S L K I K V Y D S S T I
 G S G A G L A D G A D L L G V L P P G T Y P G D F S R D T H F L H L R S A S L G S Q H L L G L P R D P S S S V S G T F G C L G G R L T I P G
 T G V S L L V P N G A I P Q G K F Y D L Y L R I N K T E S T L P L S E G S Q T V L S P S V T C G P T G L L L C R P V V L T V P H C A E V I A
 G D W I F Q L K T Q A H Q G H W E E V V T L D E E T L N T P C Y C Q L E A K S C H I L L D Q L G T Y V F T G E S Y S R S A V K R L Q L A I F
 A P A L C T S L E Y S L R V Y C L E D T P A A L K E V L E L E R T L G G Y L V E E P K T L L F K D S Y H N L X L S L H D I P H A H W R S K L
 L A K Y Q E I P F Y H V W N G S Q K A L H C T F T L E R H S L A S T E F T C K V C V R Q V E G E G Q I F Q L H T T L A E T P A G S L D A L C
 S A P G N A A T T Q L G P Y A F K I P L S I R Q K I C N S L D A P N S R G N D W R L L A Q K L S M D R Y L N Y F A T K A S P T G V I L D L W
 E A R Q Q D D G D L N S L A S A L E E M G K S E M L V A M T T D G D C


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:

SgfI-MluI

Cloning Scheme:



Plasmid Map:


ACCN: NM_022207

ORF Size: 2835 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_022207.1](#), [NP_071543.1](#)

RefSeq Size: 2838 bp

RefSeq ORF: 2838 bp

Locus ID: 60630

UniProt ID: [O08722](#)

Cytogenetics: 20q11

MW: 103.5 kDa

Gene Summary: a netrin receptor involved in neuronal migration [RGD, Feb 2006]