

Product datasheet for **RR207865L3V**

Fgf11 (NM_130816) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Fgf11 (NM_130816) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Fgf11
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_130816
ORF Size:	675 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR207865).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_130816.1 , NP_570829.1
RefSeq Size:	678 bp
RefSeq ORF:	678 bp
Locus ID:	170632
Cytogenetics:	10q24
Gene Summary:	member of the fibroblast growth factor (FGF) family; implicated in activities including embryonic development, cell growth, morphogenesis, tissue growth, and tumor growth and invasion [RGD, Feb 2006]



[View online »](#)