

Product datasheet for **RR207866L3V**

Hpca (NM_017122) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Hpca (NM_017122) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Hpca
Synonyms:	MGC105450
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_017122
ORF Size:	579 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR207866).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_017122.1 , NP_058818.1
RefSeq Size:	1561 bp
RefSeq ORF:	582 bp
Locus ID:	29177
UniProt ID:	P84076
Cytogenetics:	5q36
Gene Summary:	calcium-binding protein of the recoverin family [RGD, Feb 2006]



[View online »](#)