

Product datasheet for **RR207909L3V**

Cav2 (NM_131914) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Cav2 (NM_131914) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Cav2
Synonyms:	MGC72322
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_131914
ORF Size:	486 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR207909).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_131914.2
RefSeq Size:	2500 bp
RefSeq ORF:	489 bp
Locus ID:	363425
UniProt ID:	Q2IBC5
Cytogenetics:	4q21
Gene Summary:	a component of caveolae; involved in cholesterol transport, signal transduction, and transcytosis, may also function as a SNARE accessory protein [RGD, Feb 2006]



[View online »](#)