

Product datasheet for **RR207932**

Bend5 (NM_001108672) Rat Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Bend5 (NM_001108672) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Bend5
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin
ORF Nucleotide Sequence:	>RR207932 representing NM_001108672 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGTACGCCCTTTGTGCGGTTCTGGAGGACAACGCTCTGCTACGCGCTGCCCGTGTCTGCGTGCCGGACT
TCAGCCCCGCTCGGACTGGATTTTGACAACCAGAAGGTGTACGCGGTGTACCGGGGCCGGAGGAGCT
GGGCGTCGAGCCGAGAGCCCCCGCTACCCCTCGGGACTGGGGCGCGTGTGCTCCACAAGGCCAG
ATCCTGGCGCTGGCAGAAGACAAATCTGACCTTGAGAACAGTGTGATGCAGAAGAAAATAAAAATCCCCA
AACTTTCTCTCAGTCATGTAGAAGATGATGGGGAGGCTAAAGACTATGGGGAAGAAGATTTACAGCTTAG
ACACATCAAGAGACCGGAGGGGCGGAAGCCTAGTGAAGCGGCCACAAGAGTATCGAGGCAGTGGTGGCC
CGGCTAGAGAGACAGAACGGCCTGAGCCTGGGCCATAGCACTTGTCCAGAAGAGGTCTTTGTGGAGGCT
CAGCGGGCACAGAGGACATGGACAGTCTGGAGGATGCTGTCTGTTCCCGGGCTCTGTATGAGGAGCTGCT
GCGAACTACCAGCAGCAACAAGAAGAGATGCGCCACCTCCAGCAGGAGCTGGAGCGGACTCGGAGGCAG
CTCGTACAGCAAGCCAAGAAGCTCAAGGAATACGGGCATTGGTGTCTGAAATGAAGGAGCTCCGAGACC
TCAACAGGAGACTCCAAGATGTGCTACTCTACGCCTTGGCAGTGGTCCCGCCATTGACTTGAAAAAGT
AAAGTCAGAATGTCTCGAGCCCGAGCCGGAGTTACGGAGCACTTCAGTGAGGAAGCAAACAGTCGTCC
TATTACCCCGCTCCTGCGCCTGTGATGGACAAGTATATCCTAGACAATGGCAAGGTCCATCTGGGAAGCG
GGATTTGGGTTGACGAGGAGAAAATGGCACCAGTACAAGTTACCAAGGAGATTCTAAGTACACGAAGAA
CTTGGCAGTCAATGATCTGGGGAACAGATGTTCTGAAAAACAGAAGTGTACAGGAGTCCGCCAAAAAAA
AAGAAGGATGCAATCCCAAAGCCACCCCTCTCTCTCATAAGCTAAGCATCGTCAGAGAGTGTGTTGATG
ACAGAATAGCACAAGAAACGGTGGATGAAACGGAAATTGCACAGAGACTCTCAAAGTCAACAAGTACAT
CTGTGAGAAAAATCATGGATATCAACAAATCTTGTAAAAATGAAGAACGAAGGGAAGCAAAATACAATTTG
CAA

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >RR207932 representing NM_001108672
Red=Cloning site Green=Tags(s)


MYAFVRFLEDNVCYALPVSCVRDFSPRSRLDFDNQKVYAVYRGPEELGVEPEPSPRTPRDWGALLLHKAQ
 ILALAEKSDLENSVMQKKIKIPKLSLHVEDDGEAKDYGEEDLQLRHIKRPEGRKPSEAAHKSIEAVVA
 RLERQNGLSLGHSTCPEEVFVEASAGTEDMDSLEDVAVPRALYEELLRNYQQQEEEMRHLQQELERTRRQ
 LVQQAQKLLKEYGALVSEMKELRDLNRRLLQDVLRLRGGSPAIDLEKVKSECLEPEPELRSTFSEEANTSS
 YYPAPAPVMDKYILDNGKVHLGSGIWWDEEKWHQLQVTQGD SKYTKNLAVMIWGTDLKNRSVTGVATKK
 KKDAIPKPLSPHKLIVRECLYDRIAQETVDETEIAQRLSKVNKYICEKIMDINKSCKNEERREAKYNL
 Q

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Restriction Sites:

SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN:	NM_001108672
ORF Size:	1263 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001108672.1 , NP_001102142.1
RefSeq Size:	1834 bp
RefSeq ORF:	1266 bp
Locus ID:	362564
Cytogenetics:	5q35
MW:	48.2 kDa