

Product datasheet for **RR207935**

Selplg (NM_001013230) Rat Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Selplg (NM_001013230) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Selplg
Synonyms:	Psgl1; Selpl
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin
ORF Nucleotide Sequence:	>RR207935 representing NM_001013230 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGTTCCACACTTCCTTCTGCTGCTGACCGTCTTGGGCCCTGGCAACAGCCTTCATCTGCAGGAGCCGT
GGGGCATGAAACCATGGAAGCCCTGGGCCCTGTGCATCTCCGGGAACGGAGACAGGTGGTTGGGGATGA
TGATTTTGAATACCCGACTACAGCGTGTATGGCACAGATCCCCAGAAACGTTGATAAGTGTCCCCAAC
ATGGTGGCGGCTCACCTTGAGCTTCGGACAATGGTGGCCGTGTTAGAGCAGAGAGTTCCGCGGGGGCTG
GAATTTCTGAGACTGTCCTGCGCAGGTTGCCACCACTGGCCCTGCTGACCCGGATACAGAAAGGCCAGC
AGTTGGGATGCTGAGCACTGAATCTGTACACAGCGGAGACTAAGCCAGTGGAGATGACTACCCGGCTA
GCACCCACGGAGGCAGAGACCTCACAGCCAGCACCCATAAAGGCAGAGACCTCACAGCCAGCACCCATAA
AGGCAGAGACCTCACAGCCAGCACCCAGGGAGGCAGAGACCTCACAGCCAGCACCCACGGAGGCAGAGAC
CTCACAGCCAGCACCCACAAAGGCAGAGACCTCACAGCCAGCACCCACAGAGGCAGAGACCTCACAGCCA
GCACCCACGGAGGTAGAGACCTCACAGCCAGCACCCACAGAGGCAGAGACCTCACAGCCAGCACCCACGG
AGGCAGAGACCTCACAGCCAGCTCCACAGAGACAGAGACCCAGCTTCCCAGGAGTACAGTTGTAGA
AAGTCTATTTACAACATCTGCGGTACAGAAAGCCCTTCTACGGAACCCACCGCCATGGAGACCCCTGTCC
ACAGAGGTAATGAGTCTGCCATCTTCTTGGCCCTTCCATTGGCCACTTACCTGACACCCGGCTGAAGA
AAGGGCTGATTGTACCCCTGGGAGTTCTTGGCCCCAACCCCGCCAGGGACCTCGGACCTCATCCCGGT
GAAGCAATGTTTCTGGCCATTCTCATCTTAGCCTCCCTGGCCACCATCTTCTTGTGTGCACGGTGGT
CTGGCCGTCGCCCTCTCCCGTAAGACGCACATGTACCCAGTGGGAACTACTCCCCACAGAGATGATCT
GCATCTCATCTTGTCTCCGAGGGGGTACGGGGCCCTGTACAGCCAACGGGGCCCTGCCAAGGT
CCAGGACCAGAAGACAGAACGCAGTGGGACCCGGATGGGGATGACCTCACCTGCACAGCTTCTCCCT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >RR207935 representing NM_001013230
 Red=Cloning site Green=Tags(s)


MFPHFLLLLTVLGPNSLHLQEPWGHETMEGLGPVHLRERRQVVGGDDDFEYDPDYSVYGTDPPELISVFN
 MVAAHLELRMTMVAVLEQRVSAGAGTSETVTAQVATTGPADPDERPAVGMSTESVTQRRLSPVEMTTRL
 APTEAETSQPAPIKAETSQPAPIKAETSQPAPREAETSQPAPTEAETSQPAPTKAETSQPAPTEAETSQP
 APTEVETSQPAPTEAETSQPAPTEAETSQPASTEETTQLPRSQVVESELFSTTSAVTEAPSTEPTAMETLS
 TEANESAIFLGPSTIGHLPDTGLKKGLIVTPGSSSLAPTPPGTSDLIPVKQCLLAILILASLATIFLVCTVY
 LAVRLSRKTHMYPVRNYSPTEMICISSLLPEGGDGAPVTANGGLPKVQDQKTERSGDPDGDLLTHSFLP

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Restriction Sites:

Sgfl-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001013230

ORF Size:	1260 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001013230.1 , NP_001013248.1
RefSeq Size:	1263 bp
RefSeq ORF:	1263 bp
Locus ID:	363930
Cytogenetics:	12q16
MW:	44.3 kDa