

## Product datasheet for **RR207955L3V**

### Enpep (NM\_022251) Rat Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Enpep (NM_022251) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Enpep
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_022251
ORF Size:	2835 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR207955).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_022251.2</a> , <a href="#">NP_071587.2</a>
RefSeq Size:	3822 bp
RefSeq ORF:	2838 bp
Locus ID:	64017
UniProt ID:	<a href="#">P50123</a>
Cytogenetics:	2q42
Gene Summary:	ectoenzyme that selectively hydrolyzes acidic amino acid residues from the amino terminus of oligopeptides [RGD, Feb 2006]



[View online »](#)