

## Product datasheet for **RR207973L4V**

### **Ptk2b (NM\_017318) Rat Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	Ptk2b (NM_017318) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Ptk2b
Synonyms:	CAKB; CAK beta; CAKbeta; Pyk2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_017318
ORF Size:	3027 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR207973).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_017318.2</a> , <a href="#">NP_059014.2</a>
RefSeq Size:	3907 bp
RefSeq ORF:	3030 bp
Locus ID:	50646
UniProt ID:	<a href="#">P70600</a>
Cytogenetics:	15p12
Gene Summary:	may phosphorylate focal adhesion kinase FAK; may play a role in actin cytoskeleton reorganization; may mediate thrombin induced astrocyte proliferation [RGD, Feb 2006]



[View online »](#)