

Product datasheet for **RR207976L4V**

Dync1li2 (NM_031026) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Dync1li2 (NM_031026) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Dync1li2
Synonyms:	Dncli2; Lic2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_031026
ORF Size:	1491 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR207976).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_031026.1 , NP_112288.1
RefSeq Size:	4300 bp
RefSeq ORF:	1494 bp
Locus ID:	81655
UniProt ID:	Q62698
Cytogenetics:	19p14
Gene Summary:	light intermediate form of cytoplasmic dynein involved in retrograde organelle transport and mitosis; may regulate the interactions between dynein and p150-glued [RGD, Feb 2006]



[View online »](#)