

## Product datasheet for **RR207986**

### Rab40b (NM\_001107076) Rat Tagged ORF Clone

#### Product data:

|                          |  |
|--------------------------|--|
| Product Type:            | Expression Plasmids  |
| Product Name:            | Rab40b (NM_001107076) Rat Tagged ORF Clone                                     |
| Tag:                     | Myc-DDK  |
| Symbol:                  | Rab40b   |
| Vector:                  | pCMV6-Entry (PS100001)   |
| E. coli Selection:       | Kanamycin (25 ug/mL)   |
| Cell Selection:          | Neomycin   |
| ORF Nucleotide Sequence: | >RR207986 representing NM_001107076<br>Red=Cloning site Blue=ORF Green=Tags(s) |

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGATGAGCTCCTTGGGCAGCCGGTGC GGCCCTACGATTTTTACTTAAGTTCCTGTTGGTGGCGACA  
GCGACGTGGCAAGGGAGAGATCCTGGCCAGCTTGCAGGACGGAGCAGCCGAATCTCCTTACGGCCACC  
GGCGGGAATTGATCACAAGACGACCACCATCTGCTAGATGGCCGAAGGGTGAAACTGCAGCTTTGGGAC  
ACCTCCGGCAAGGGAGATTTGTACCATTCCGCTCGTACTCCAGAGGTGCACAGGGTGTGATCCTGG  
TCTATGACATTGCCAACCGATGGTCCTTTGATGGCATTAAATCGGTGGATTAAGGAGATTGATGAGCATGC  
CCCTGGGGTTCCTAAAATCTTGGTGGGAAACCGCTACACCTGGCATTCAAGCGGCAGGTACCCACCGAG  
CAGGCCCAGGCCCTACGCTGAGCGCCTGGGTGTACCTTCTTGAAGTTAGCCCTCTGTGCAACTTCAACA  
TTACAGAGTCCTTACGGAGCTGGCAAGGATAGTGTCTGCTGCGGCATGGGATGGACCGTCTCTGGAGGCC  
AAGCAAAGTACTGAGTTTGAAGACCTCTGCTGCCGTGCTGTGGTATCCTGCACACCTGTGCACCTGGTA  
GACAAGCTCCCTCTCCCTGTTGCCTTAAGAAGCCACCTTAAGTCCTTCTCTATGGCCAACGGCCTGAATA  
CCAGGATGATGCATGGCCGCTCCTACTCCCTCACTGCCAACTCCAGCCATAAGAGGAACAGCTTCCGAAA  
AGTCAGGACCATCCGCCACCCAGAGTCCACCCAGAACTGTGCCAGGAACAGCTGCAAAATTTCT

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)


**Protein Sequence:** >RR207986 representing NM\_001107076  
 Red=Cloning site Green=Tags(s)

MMSLGLSPVRAYDFLLKFLLVGDSVDVGKGEILASLQDGAESPYPAGIDHKTTTILLDGRRVKLQLWD  
 TSGQGRFCTIFRSYSRGAQGVILVYDIANRWSFDGINRWIKEIDEHAPGVPKILVGNRHLAFKRQVPT  
 QAQAYAERLGVTFEVSPLCNFNITESFTELARIVLLRHGMDRLWRPSKVLSLQDLCCRAVVSCTPVHLV  
 DKLLPLVALRSHLKSFSMANGLNTRMMHGRSYSLTANSSHKRNSFRKVRTIRPPQSPPRNCARNSCKIS


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**



**Plasmid Map:**



**ACCN:** NM\_001107076

**ORF Size:** 837 bp

|                               |   |
|-------------------------------|---|
| <b>OTI Disclaimer:</b>        | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>  |
| <b>OTI Annotation:</b>        | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.  |
| <b>Components:</b>            | The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).  |
| <b>Reconstitution Method:</b> | <ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol> |
| <b>RefSeq:</b>                | <a href="#">NM_001107076.1</a> , <a href="#">NP_001100546.1</a>   |
| <b>RefSeq Size:</b>           | 1861 bp   |
| <b>RefSeq ORF:</b>            | 840 bp  |
| <b>Locus ID:</b>              | 303754  |
| <b>Cytogenetics:</b>          | 10q32.3   |
| <b>MW:</b>                    | 31.3 kDa  |