

Product datasheet for **RR208005L3V**

Dcn (NM_024129) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Dcn (NM_024129) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Dcn
Synonyms:	MGC94682
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_024129
ORF Size:	1062 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR208005).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_024129.1 , NP_077043.1
RefSeq Size:	1639 bp
RefSeq ORF:	1065 bp
Locus ID:	29139
UniProt ID:	Q01129
Cytogenetics:	7q13
Gene Summary:	binds the N-terminal region of collagen VI; may play a role in extracellular matrix assembly [RGD, Feb 2006]



[View online »](#)