

## Product datasheet for **RR208018**

### Zmat2 (NM\_001135582) Rat Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Zmat2 (NM\_001135582) Rat Tagged ORF Clone  
**Tag:** Myc-DDK  
**Symbol:** Zmat2  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**ORF Nucleotide Sequence:** >RR208018 representing NM\_001135582  
**Red=Cloning site Blue=ORF Green=Tags(s)**

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGGCGTCTGGCAGCGGGACAAAACTTGGACTTTCGCCGAAAGTGGGACAAAGATGAATATGAGAAGC  
TCGCGGAGAAGAGACTCACAGAAGAGAGAGAAAAGAAGGATGGAAAACAGTGCAGCCAGTTAAGCGGGA  
GCTTCTACGGCACAGGGATTACAAGGTGGACTTGGAGTCCAAGCTTGGGAAGACAATTGTTATACCAAG  
ACCACCCACAGTCTGAAATGGGAGGCTACTACTGCAATGTCTGTGACTGTGTGGTGAAGGACTCCATTA  
ACTTCTGGATCACATTAATGGAAGAACAATCAGCGAAACCTGGGTATGTCTATGCGCGTGAACGTTT  
CACTCTGGATCAGGTAAGAAACGTTTTGAGGTTAACAAAGAAGAAGATGGAGGAGAAGCAGAAAGATTAT  
GATTTTGAGGAAAGGATGAAGGAGCTCAGAGAAGAGGAGGAAAAGGCTAAAGCATAACAAGAAAGAGAAAC  
AGAAGGAGAAGAAAGGAGGGCCGAGGAGGACTTGACTTTTGAGGAGGAGGATGAAATGGCAGCTGTGAT  
GGGCTTCTGGCTTTGGTTCTACCAAGAAGATTAC

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >RR208018 representing NM\_001135582  
**Red=Cloning site Green=Tags(s)**

MASGSGTKNLDFRRKWDKDEYEKLAEKRLTEEREKKDGKPVQPVKRELLRHRDYKVDLESKLGKTIVITK  
TTPQSEMGGYYCNVDCVVKDSINFLDHINGKKHQENLGMMSRVERSTLDQVKKRFEVNNKKMEEKQKDY  
DFEERMKELREEEKAKAYKKEKQKEKRRRAEEDLTFEEDEMAAVMGFSGFGSTKKS

**TR**TRPLEQKLISEEDLAANDILDYKDDDDKV

**Restriction Sites:** SgfI-MluI



**Cloning Scheme:**


**ACCN:** NM\_001135582

**ORF Size:** 597 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
  2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
  3. Close the tube and incubate for 10 minutes at room temperature.
  4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
  5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001135582.1](#), [NP\\_001129054.1](#)

**RefSeq Size:** 1809 bp


**RefSeq ORF:** 600 bp

**Locus ID:** 679898

**Cytogenetics:** 18p11

**MW:** 23.6 kDa

Product images:



Circular map for RR208018