

Product datasheet for **RR208020**

Pdzk1 (NM_031712) Rat Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Pdzk1 (NM_031712) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Pdzk1
Synonyms:	Clamp
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

ORF Nucleotide Sequence:

>RR208020 representing NM_031712
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGCC**

ATGGCCTCCACCTTCAACCCAGAGAGTGTAAATTGTCCAAAAAGAGGGGCAGAACTATGGCTTCTTCC
 TCCGATTGAAAAGGACACTGACGGTCACTGGTCCGGTGATTGAGGAGGGAGCCCGCAGAGAAGGC
 AGGGCTCTAGACGGCGACAGAGTGCTCAGGATCAACGGGGTCTTTGTTGACAAGGAAGACATGCACAG
 GTGGTGGATCTGGTCAGAAAAGAGCGGAATTCAGTGACCTGCTGGTTCTGGATGGAGACTCCTATGAGA
 AGGCTGTAAAACATCAGGTGGACTTGAAGAAGTGGATCAAAGCCCAAGGGAGCCAGCTCTGAATGAAAA
 GAAACCCGACCTTGGGATGAATGGAGGGGTGGAGACGTGTGCCACGCCAGGCTCTGCTACCTGGTGAAG
 GAGGGCAACAGCTTTGGCTTCTCTGAAAATCAAGGTAAAAGGGCGTGTCTGACTGATATAA
 CACCTCAGGGCGTGGCCATGAAAGCTGGTGTCTGGCTGATGATCACTTGATTGAAGTGAATGGAGAAAA
 TGTAGAAAACCCAGTCATGAGGAAGTGGTGGAAAAGGTGACAAAGTCAGGAAGCCGTATCATGTTCCCTC
 CTTGTGGACAAAGAACTGCCAGGTGCCATAGTGAGCAGAAAACACCGTTCAAGAGGGAGACAGCCAGTT
 TGAAACTGCTGCCCCACAGCCCCGGTAGTAGTGATCAAGAAGGGCAGCAATGGCTACGGTTTCTATCT
 GAGGGCGGCCCTGAACAGAAAGGTCAAATCATCAAGGACATAGAACCTGGAAGCCAGCAGAGGCAGCT
 GGCTTGAAGAACAATGACTTGGTAGTTGCTGTCAATGGCAGTCTGTGGAAGCTCTTGATCATGACAGTG
 TGGTGGAAATGATTAGAAAACGGTGGAGACCAGACTACTCTGTTGGTGTGGACAAGGAGGAGCAGAGAA
 CTATAGCCTGGCTCGTTTCTCTCCACTTCTCTACTGCCAAAGTCAAGAACTGCCTAATGGTTCTGTGAAG
 GAAGCCCCAGCTCCCATCTCTGCTCCTCTGGAGGCCAGGCTCAGCTACCACAGAGGACGTGGGGGATC
 ATAAACCAAGCTCTGCAGGCTGATTAAGAGGACGACAGCTACGGCTTTCACCTGAACGCCATTCGGGG
 CCAGCCAGGCTCCTTTGTCAAGGAGGTACAGCAGGGTGGCCCTGCTGACAAGGCAGGGTGGAGAATGAG
 GACATCATCATTGAAGTGAATGGGGAGAATGTGCAAGATGAACCCATGACAGAGTGGTGGAGAGATCA
 AGAGCAGCGGGGAGCACGTCACTCTGTTGGTCTGTGAAAAGTGGCTACAGCTATTTTCAAGCTAAGAA
 AATCCCCATCCTTTCTTCTGGCTGACCCCTGGTGGCTGGTCTGATGAAAAAGGAGAGACTGAGCAT
 GACTCAGCAGAGTCCACGAAAGACTCTCCACCCAGCAAGAGACCGAACTCTCAGCGCAGCCTCCCACT
 CTTCTTCCAACCTGAAGACACGGTGATG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:

>RR208020 representing NM_031712
 Red=Cloning site Green=Tags(s)


MASTFNPRECKLSKKEGQNYGFFLRIEKDTDGHVLRVIEEGSPAELLDGDRVLRINGVFDKEEHAQ
 VVDLVRKSGNSVTLVLDGDSYEKAVKHQVDLKDQSPREPALNEKKPDLGMNGGVETCAQPRLCYLK
 EGNSFGFSLKTIQGGKGVFLDITPQGVAMKAGVLADDHLIEVNGENVENASHEEVVEKVTKSGSRIMFL
 LVDKETARCHSEQKTPFKRETASLKLPHQPRVVVIKKGSNGYGFYLRAGPEQKQI IKDIEPGSPAEAA
 GLKNNDLVAVNGESVEALDHSVEMIRNGDQTLLVLDKEADRIYSLARFSPLLYCQSQELPNGSVK
 EAPAPISAPLEAPGSATTEDVGDHDKPKLRLIKEDDSYGFHLNARIRGPGSFVKEVQGGPADKAGLENE
 DIIIEVNGENVQDEPYDRVVERIKSSGEHVTLVCGKVAYSFYQAKKIPILSSLADPLVAGPDEKGETEH
 DSAESTKDSSHPARDRTL SAASHSSNSSEDTVM

TRTRPLEQKLI SEEDLAANDILDYKDDDDKV


Restriction Sites:

Sgfl-MluI

Cloning Scheme:



Plasmid Map:



ACCN:	NM_031712
ORF Size:	1569 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_031712.1 , NP_113900.1
RefSeq Size:	2005 bp
RefSeq ORF:	1572 bp
Locus ID:	65144
UniProt ID:	Q9JJ40
Cytogenetics:	2q34
MW:	56.9 kDa
Gene Summary:	binds the N-terminal PDZ domain of scavenger receptor class B type I; may play a role in intracellular transport of cholesterol esters [RGD, Feb 2006]