

Product datasheet for **RR208020L3V**

Pdzk1 (NM_031712) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Pdzk1 (NM_031712) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Pdzk1
Synonyms:	Clamp
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_031712
ORF Size:	1569 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR208020).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_031712.1 , NP_113900.1
RefSeq Size:	2005 bp
RefSeq ORF:	1572 bp
Locus ID:	65144
UniProt ID:	Q9J140
Cytogenetics:	2q34
Gene Summary:	binds the N-terminal PDZ domain of scavenger receptor class B type I; may play a role in intracellular transport of cholesterol esters [RGD, Feb 2006]



[View online »](#)