

Product datasheet for **RR208022**

Foxs1 (NM_001012091) Rat Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Foxs1 (NM_001012091) Rat Tagged ORF Clone
Tag: Myc-DDK
Symbol: Foxs1
Synonyms: Fkhl18
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RR208022 representing NM_001012091
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGCAGCAGCCACCCACAGCAGAATCCTTGCTCCTTCTACAGAGCCCAGCAAGCCCCCTTACAGCTACA
TAGCCCTGATTGCCATGGCTATCCAGAGTTCACCGGGTCAGCGCGCCACCCTCAGTGGCATCTACCCTG
CATCATGGGCGCTTTGCATTTTACCGGCACAACCGGCCAGGCTGGCAAACAGCATCCGCCATAACCTG
TCTCTCAACGAGTGCTTCGTCAAGGTGCCTCGTGATGACCGCAAGCCAGGCAAGGGCAGCTACTGGACCC
TGGATCCAGACTGCCACGACATGTTCCAGCACGGCAGCTTCTCCGTCCGCCGAGACGCTTACCAAGCG
CACAGGTGCCAGGGCACCAAGGGCCCTGTCAAGGCAGACCACAGACCCCTTAGAGCCACCAGCCCGGAT
CAAGGAGCCCCAACACCACAACCTGGCAGACTGTGCCATTCCCACCAGAGGTGCCAAACCCCAAGGGCT
TTGGGGTCTGATGGGGTCATTGCCAGCCAACATGTGCCAACCAACAGTGACACCAGGCCCCAGCTGCC
CACTGGACCCAAAGACATGTGCTCGGCAAAGTCTGGAGGTCCAAGGGAGCTTTCTGAAGCCACCTCGCCC
TCACCATGCCCGCATTCGGCTTCTCATCTGCCTTCTCTGAAGCTGAGAGTCTCGGCAAGGCCCCACAC
CCAGTGTGGCTCCGGAGAGCATCGGGAGCAGTACCAGTCCGGATGCAGACGCTGAATTTCTGTATGGG
GGCTGACCCAGGCCTTGAACACCTTTGACCTCTGCAGTCCGCACTCCTGGATCCTCAACACCTTCGGCC
TCTCACCGGGCTCCACTGCCCTGCCAGTGACTCCAAGGAACCGTGGGTGCCTGGGGGCTTCCCTGTCC
AGGGAGTTCTGGCTATCCACTGGGATTACCCCTTGTCTTTACCGGACACCGGGAATGTTCTTCTCGA
G

ACGCGTACGCGGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)


Protein Sequence: >RR208022 representing NM_001012091
 Red=Cloning site Green=Tags(s)

MQQPPTAESLAPSTEPSKPPYSYIALIAMAIQSSPGQRATLSGIYRYIMGRFAFYRHNRPGWQNSIRHNL
 SLNECFVKVPRDDRKPGKGSYWTLPDCHDMFQHGSLRRRRRF TKRTGAQGTGKPVKADHRPLRATSPD
 QGAPNTTTGRLCPFPPEVPNPKFGGLMGSLPANMCPTTSDTRPQLPTGPKDMCSAKSGGPRELSEATSP
 SPCPAFGFSSAFSEAESLGKAPTPSVAPESIGSSYQCRMQTLNFCMGADPGLHLLTSAVATPGSSTPSA
 SHRAPLPLPADSKEPWPVGGFPVQGGSGYPLGLPPCLYRTPGMFFFE


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001012091

ORF Size: 981 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001012091.1 , NP_001012091.1
RefSeq Size:	1300 bp
RefSeq ORF:	984 bp
Locus ID:	311547
Cytogenetics:	3q41
MW:	35.2 kDa