

Product datasheet for **RR208055L4V**

Apeh (NM_012500) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Apeh (NM_012500) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Apeh
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_012500
ORF Size:	2196 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR208055).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_012500.1 , NP_036632.1
RefSeq Size:	2381 bp
RefSeq ORF:	2199 bp
Locus ID:	24206
UniProt ID:	P13676
Cytogenetics:	8q32
Gene Summary:	catalyzes the removal of N alpha-acetylated amino acid residues from peptides [RGD, Feb 2006]


[View online »](#)