

Product datasheet for **RR208056L3V**

Fetub (NM_053348) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Fetub (NM_053348) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Fetub
Synonyms:	Fet; IRL685; Pp63
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_053348
ORF Size:	1134 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR208056).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_053348.2 , NP_445800.2
RefSeq Size:	1423 bp
RefSeq ORF:	1137 bp
Locus ID:	83928
UniProt ID:	Q9QX79
Cytogenetics:	11q23
Gene Summary:	inhibits insulin receptor tyrosine kinase activity, protease inhibitory activities; developmentally regulates calcium metabolism and osteogenesis [RGD, Feb 2006]



[View online »](#)