


Product datasheet for **RR208071L3**

Adnp (NM_022681) Rat Tagged Lenti ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Adnp (NM_022681) Rat Tagged Lenti ORF Clone
Tag:	Myc-DDK
Symbol:	Adnp
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR208071).
Restriction Sites:	Sgfl-MluI
Cloning Scheme:	

Cloning sites used for ORF Shuttling:



ACCN:	NM_022681
ORF Size:	3309 bp



View online »

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_022681.2 , NP_073172.2
RefSeq Size:	5608 bp
RefSeq ORF:	3312 bp
Locus ID:	64622
Cytogenetics:	3q42
Gene Summary:	This gene encodes a member of a protein family characterized by nine zinc finger motifs followed by a homeobox domain. In vitro studies demonstrate that the encoded protein interacts with the brahma-related gene1-associated or hBRM factors (BAF) gene expression regulating complex, components of the protein translation machinery, and microtubule-associated proteins. This gene has been implicated in neuroprotection through various processes that include chromatin remodeling, splicing, cytoskeletal reorganization, and autophagy. Homozygous mutant knockout mice display embryonic lethality with defects in neural tube closure. [provided by RefSeq, Nov 2016]