

Product datasheet for **RR208118**

Snap47 (NM_199389) Rat Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
 Product Name: Snap47 (NM_199389) Rat Tagged ORF Clone
 Tag: Myc-DDK
 Symbol: Snap47
 Synonyms: RGD735194; SNAP-47
 Vector: pCMV6-Entry (PS100001)
 E. coli Selection: Kanamycin (25 ug/mL)
 Cell Selection: Neomycin
 ORF Nucleotide Sequence: >RR208118 representing NM_199389
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGAGCAGTGATGTACGAGTTCACACCTGGCCTTGTTCTACTATCTGGACCTTGAGAAGCAATGGCTTC
 CTGGAAAACCTGACCTTGACACCTCGTTCCTGACTGAAGTTCAGTGTGGACAAAGCTGAAGAGGCTCTTGTGG
 CCTGCCCTTTCTAGCATCACTGAGATCAGGAAAGAGACATCGCTTTTATCTTCAGTGGATCACAGTC
 CTGGAGAAGGGCCAGGCCAAGCACTGGTTCAGCTCACTGCGGCCTAGTCGAAATGTGGTCTTCAACATCA
 TCGAGCATTCTGGAGGGAGCTGCTACTGTCCCAGCCAGGCACTGCTGCCAATCCAACCTCCCCATGAC
 CAGGGGGCAGGAGCTGATGGGACTGATGGCGAGTTCAGAGACGCATGGAAGACACAGCAAAGGTCCTT
 CACCACCAGGGCGAGCAGCTAGACAGTGTATGAGGGCCTGGAGAAGATGGAGTCCGATCTAGATGTGG
 CTGACAGATTGCTGACAGAGTTGGAGACTCCTTCTGTTGGCCTTTTGTAGTCTAAGTTTTGGAAGACGCC
 AGTGAAGCCAAGCCAGGGACAGTGTCTGTAGCCCTTGTGAGCCCTTTGGAAGGAAGGGTTGTG
 ATCAGGGTCCCTGCTGTTGTTTCTCAGAGAACAGAGTCTCATTCCAAGCCAGGAAAACCTACTGTCCTGG
 TCTCCGGTTTGGAAATCCATGACTCCAGTTCATTGCTCCTGCACAGATTTGAGAAGGAAGATGTTGATGA
 TATCAAGTTTCACTCTCCTTATGAGGTACGATACGCCAGCGGTTTCATCGGCAAGCCAGAGCTGGCTTAC
 CGGCTGATATCTGCCAAGATGCCAGAGGTCACTCCATTTTAGAAGTGCAGTTTCAGCAAGAAGATTGAGC
 TCCTGGAAGACGCCCTGGTGTCTCAGGAGCAGGGCCGTGCTTCGCCTGCCGAGGGCGGCTGCTCCATTTCG
 GCATGCAGCATCCAGACTGATGGGCTGCACTACACACCAAGAGCTGCCACAGGAGGACAGGAAGGCCAG
 CAGTCCGAGCTGCAGAAGGACTGGCCACTACTCTGAGGGGGAAGCCAGGAGCTGACACAGATCTGA
 GCAAACGAAGGGCCTCGCACTGGACACCGAGGACAGCTGGAACGCCAGGATGCAGCTCTGGATGGCAT
 CACCGTAGCTGTGGACCGGGCAACCCTGACTGTGGACAAGCATAACCGTCGTATGAGAAAGCTGTTG

ACGCGTACGCGGCCGCTCGAGCAGAAAACCTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >RR208118 representing NM_199389
 Red=Cloning site Green=Tags(s)


MSSDVRVHTWPCSYLLDLEKQWLPGLTLTPRSLKFSVDKAEVLVGLPLSSSITEIRKETSLFIFSAITV
 LEKGQAKHWFSSLRPSRNVVFNIEHFWRELLLSQPGTAAAPTSPMTRGQELMGLMASSQRRMEDTAKVL
 HHQGEQLDSVMRGLKEMESDLDVADRLLTELETPSWWPFSAKFWKTPVEAKPRDSVSVAPCEPFGKEGVV
 IRVPAVVSQRTESHKPGKLTVLVSGLEIHDSSLLLLHRFEKEDVDDIKVHSPYEVSIQRFIGKPDVAY
 RLISAKMPEVPIILEVQFSKKIELLEDALVLRSGRASPAGGCSIRHAASRLMGCTTHQELPTGGQEQQ
 QSQLQKDWPLLEGEAQELTQILSKLKGLALDTEAELERQDAALDGI TVAVDRATLTVDKHNRMRKLL

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Restriction Sites:

SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN:

NM_199389

ORF Size:	1257 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_199389.1 , NP_955421.1
RefSeq Size:	1822 bp
RefSeq ORF:	1260 bp
Locus ID:	303183
UniProt ID:	Q6P6S0
Cytogenetics:	10q22
MW:	47 kDa
Gene Summary:	May play a role in intracellular membrane fusion.[UniProtKB/Swiss-Prot Function]