

Product datasheet for **RR208121**

Cdk2ap2 (NM_001109498) Rat Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Cdk2ap2 (NM_001109498) Rat Tagged ORF Clone
Tag: Myc-DDK
Symbol: Cdk2ap2
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RR208121 representing NM_001109498
Red=Cloning site **Blue**=ORF **Green**=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGTCCTACAAACCTATCGCGCCCGCCCCAGCAGTACCCCGGGCTCCAGCACCCCGGGCCAGGCACCC
CGGTCCCTACAGCCGGAAGTGTCCATCGCCATCGGGCTCAGTGCCAGGAGCTGCAGCTCCTTTTCAGACC
ACTGTTTAACGACTTTGGACCGCCGTCATGGGGTATGTACAGGCAATGAAGCCACCCGGTTCCAGGGC
TCTCAGAGCACCTACACGACCTGCTGTCTGTATAGAGGAGATGGGCAAAGAGATCCGACCCACCTATG
CTGGTAGCAAGAGTGCCATGGAGCGCCTGAAGAGAGGTATCATCCATGCTCGGGCACTAGTCAGAGAGTG
CTTGCCAGAGACAGAACGCAATGCCCGCAGC

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RR208121 representing NM_001109498
Red=Cloning site **Green**=Tags(s)

MSYKPIAPAPSSTPGSSTPGPGTPVPTAGSVSPSGSVPGAAAPFRPLFNDFGPPSMGYVQAMKPPGSQG
SQSTYTDLLSVIEEMGKEIRPTYAGSKSAMERLKRGIHARALVRECLAETERNART

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI



[View online »](#)

OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001109498.1</u> , <u>NP_001102968.1</u>
RefSeq Size:	907 bp
RefSeq ORF:	384 bp
Locus ID:	688405
Cytogenetics:	1q43
MW:	13.2 kDa