

Product datasheet for **RR208146L4V**

Cct3 (NM_199091) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Cct3 (NM_199091) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Cct3
Synonyms:	MGC72950
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_199091
ORF Size:	1635 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR208146).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_199091.2 , NP_954522.1
RefSeq Size:	1981 bp
RefSeq ORF:	1638 bp
Locus ID:	295230
UniProt ID:	Q6P502
Cytogenetics:	2q34
Gene Summary:	one of eight subunits of CCT, an ATP-ase and chaperonin protein which aids in the ATP-dependent folding of actin and tubulin in the cytosol [RGD, Feb 2006]



[View online »](#)