

Product datasheet for **RR208247L4V**

Bckdha (NM_012782) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Bckdha (NM_012782) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Bckdha
Synonyms:	BCKDA; E1a
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_012782
ORF Size:	1338 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR208247).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_012782.1 , NP_036914.1
RefSeq Size:	1690 bp
RefSeq ORF:	1341 bp
Locus ID:	25244
UniProt ID:	P11960
Cytogenetics:	1q21
Gene Summary:	plays a role in the conversion of 3-methyl-2-oxobutanoate and lipoamide to S-(2-methylpropanoyl)dihydrolipoamide and CO2 [RGD, Feb 2006]



[View online »](#)