

Product datasheet for **RR208260L4V**

Adcy2 (NM_031007) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Adcy2 (NM_031007) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Adcy2
Synonyms:	AC2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_031007
ORF Size:	3270 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR208260).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_031007.1 , NP_112269.1
RefSeq Size:	4008 bp
RefSeq ORF:	3273 bp
Locus ID:	81636
UniProt ID:	P26769
Cytogenetics:	1p11
Gene Summary:	enzyme that catalyzes the formation of the secondary messenger cyclic adenosine monophosphate [RGD, Feb 2006]



[View online »](#)