

Product datasheet for **RR208324L3V**

Uchl3 (NM_001110165) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Uchl3 (NM_001110165) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Uchl3
Synonyms:	RGD1561196
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001110165
ORF Size:	690 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR208324).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001110165.1 , NP_001103635.1
RefSeq Size:	914 bp
RefSeq ORF:	693 bp
Locus ID:	498560
UniProt ID:	Q91Y78
Cytogenetics:	15q21
Gene Summary:	mouse homolog is an enzyme involved in proteolytic processing of polymeric ubiquitin [RGD, Feb 2006]



[View online »](#)