

Product datasheet for **RR208330L3V**

Ccl3 (NM_013025) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Ccl3 (NM_013025) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Ccl3
Synonyms:	MIP-1a; Scya3
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_013025
ORF Size:	276 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR208330).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_013025.2 , NP_037157.2
RefSeq Size:	755 bp
RefSeq ORF:	279 bp
Locus ID:	25542
UniProt ID:	P50229
Cytogenetics:	10q26
Gene Summary:	mediates monocyte and neutrophil chemotaxis; may play a role in the pathogenesis of acute lung injury [RGD, Feb 2006]



[View online »](#)