

Product datasheet for **RR208344**

Lsp1 (NM_001025420) Rat Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Lsp1 (NM_001025420) Rat Tagged ORF Clone
Tag: Myc-DDK
Symbol: Lsp1
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RR208344 representing NM_001025420
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGCAGAGGCTGCCAGCGATCCCAGATGCCAGGAGCAGGAAGAGCTAGATGCAGAAGACGACAAAGGGC
TCACAACCTCAGTGGAGAGAAGAAGATGAAGAAGAGGCCGCCCGTGAGCAGCGCCAACGGGAGCGAGAGAG
GCAGCTGCAGACCAGGACAAAGATGAAGAAGATGAAGGTGGCCATTTCCCTGGAACAGCCAGGGCAGCAG
GCACTCGTTAGCCTGAAGTCCTCTGAACTGGATGAAGATGAGGGTTTTGGTGACTGGTCTCAAAGCAGG
AGCCACCAAGGCAGCAATCCTGGGAGAATGAGGGGACTGCAGAGGGTCCGAACCCCTCTCCAAGTGAGAG
TGCAGAGGAGAAACAAACAGAGGACAGTTCTACCAAGCTAAAGTCCACTTGGAGGAGTCAAACCTGAGC
CACAGTGACCCCAATATAGAGGATGATGTCGGGGTCTGGGGAGACAGAAGAGCATCTGATAAGTCATC
AGGCCAGGACCCCGCCCTTTGGCCTTGAAGACTGCTGAGCTGAGCTCACCTCCCCTGAGCCCTAC
CGTCAAACCTGGCTGATAGGACAGAATCCCTGAACCGCTCCATTCAAAGAGCAACAGTGTGAAGAAGTCT
CAACCAACCTTGCCATTTCCACCATCGATGAGCGCTGCAGCAGTACACCCAGGCCACTGAGTCTGCTG
GCCGAACCTCCAAGCTGTCCCGCCAGCCCTCCATAGAGCTGCCAGCATGGCTGTAGCCAGTACCAAGAC
TCTTTGGGAGACAGGAGAGGTGCAGACTCAGTCTGCTTCCAAGACACCCTCCTGCCAGGATATAGTAGT
GGAGACATGAGCAAGAAAAGTCTGTGGGAACAGAAAGGAGGCTCCAAGATCTCATCCACCATCAAGAGCA
CCCCATCTGGAAGAGGTACAAGTTTGTGGCCACTGGACACGGGAAGTACGAGAAAGTGCTCGTGGATGA
GGGCTCAGACCA

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >RR208344 representing NM_001025420
 Red=Cloning site Green=Tags(s)

MAEAAADPRCQEELDAEDDKGLTTQWREDEEEAAREQRQREERERQLQDQDKDEEDEGGHFLEQPGQQ
 ALVSLKSSSELDDEDEFGDWSQKQEPQRQSWENEGTAESESPSESAEEKQTEDSSHQAKVHLEESNLS
 HSDPNIEDDVGGSGETEELISHQARTPSPLALEDTAELSSPPLSPTVKLADRTESLNRSIQKSNVKKKS
 QPTLPISIDERLQQYTQATESAGRTPKLSRQPSIELPSMAVASTKTLWETGEVQTQSASKTPSCQDIVA
 GDMSKSLWEQKGGSKISSTIKSTPSGKRYKFVATGHGKYEKVLVDEGSAP


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001025420

ORF Size:	993 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001025420.1 , NP_001020591.1
RefSeq Size:	1257 bp
RefSeq ORF:	996 bp
Locus ID:	361680
Cytogenetics:	1q41
MW:	36.6 kDa