

Product datasheet for **RR208360L4V**

Tesk1 (NM_031578) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Tesk1 (NM_031578) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Tesk1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_031578
ORF Size:	1884 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR208360).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_031578.1 , NP_113766.1
RefSeq Size:	3581 bp
RefSeq ORF:	1887 bp
Locus ID:	29460
UniProt ID:	Q63572
Cytogenetics:	5q22
Gene Summary:	serine/threonine protein kinase; phosphorylates myelin basic protein (MBP) and histones in vitro; may play a central role in spermatogenesis [RGD, Feb 2006]



[View online »](#)