

Product datasheet for **RR208362L4V**

Msi1 (NM_148890) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Msi1 (NM_148890) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Msi1
Synonyms:	Msi1h
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_148890
ORF Size:	1086 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR208362).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_148890.1 , NP_683688.1
RefSeq Size:	1089 bp
RefSeq ORF:	1089 bp
Locus ID:	259272
UniProt ID:	Q8K3P4
Cytogenetics:	12q16
Gene Summary:	RNA binding protein; displays increased expression in reactive astrocytes and in the subgranular zone (SGZ) of the hippocampal dentate gyrus after ischemia [RGD, Feb 2006]



[View online »](#)