

Product datasheet for **RR208369L3V**

Igfbp5 (NM_012817) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Igfbp5 (NM_012817) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Igfbp5
Synonyms:	IGF-BP5
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_012817
ORF Size:	813 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR208369).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_012817.1 , NP_036949.1
RefSeq Size:	1630 bp
RefSeq ORF:	816 bp
Locus ID:	25285
UniProt ID:	P24594
Cytogenetics:	9q33
Gene Summary:	may be involved in ovarian folliculogenesis and may play a modulatory role in type I fiber-dominated muscles [RGD, Feb 2006]



[View online »](#)