

## Product datasheet for **RR208375L4V**

### Akt3 (NM\_031575) Rat Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Akt3 (NM_031575) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Akt3
Synonyms:	Pkbg
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_031575
ORF Size:	1362 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR208375).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_031575.1</a> , <a href="#">NP_113763.1</a>
RefSeq Size:	1548 bp
RefSeq ORF:	1365 bp
Locus ID:	29414
UniProt ID:	<a href="#">Q63484</a>
Cytogenetics:	13q25
Gene Summary:	kinase that associates with both protein kinase C subspecies and beta gamma subunits of G proteins [RGD, Feb 2006]



[View online »](#)